



# Leica EM PACT2

High Pressure Freezer


# Leica EM RTS

Rapid Transfer System

# High Pressure Freezing for Everybody

The Leica EM PACT2 high pressure freezer serves the needs of molecular and cell biologists and all researchers who want an “in vivo” impression of their cellular structures and functions in question – without the artefacts of chemical fixation but with the high resolution information of EM immunocytochemistry, frozen hydrated sections and freeze fracturing.

The Rapid Transfer System EM RTS allows correlative LM/EM experiments, taking a live specimen from a light microscope (e.g. a confocal microscope) to freezing in less than 5 seconds. In the same way, time resolved experiments are possible. Safety and reproducibility for the specimen are increased while operator mistakes are reduced.



Structural details of the *C. elegans* head in cross-section.  
T. Müller-Reichert (MPI-CBG, Dresden, Germany) and Kent McDonald (University of California, Berkeley, USA)





Leica Design by W. Hölbl

# Leica EM PACT2 – High Pressure Freezer



## Perfect results

- High cooling rates by strong jet of LN<sub>2</sub>
- 120 samples can be frozen per hour
- Variety of specimen carriers available for all purposes

## Easy to install


- Compact and mobile
- Standard electrical supply
- Compressor on EM PACT2 trolley

## Easy to use

- Touch-sensitive color screen, menu prompts
- Bayonet loading device for automatic orientation
- Specimen ejection into LN<sub>2</sub> bath
- LED bath illumination
- 600 ms temperature/pressure curve displayed for each run
- Internal memory for 8 000 freezing runs
- Data download on memory stick

## Safe and convenient

- Low LN<sub>2</sub> consumption
- Dewar with drain outlet
- Universal hydraulic system – can even use water
- Low noise
- Maintenance free “long life pods”



### Function: All You Need is LN<sub>2</sub>...

The specimen is locked (1) and set under pressure (2) just before freezing (3) via the synchronization mechanism (4). Only the tissue in the specimen holder is under pressure and not the complete specimen chamber. As any fluid can be used in the universal hydraulic system, toxic vapors (e.g. ethanol) are avoided.

# Leica EM RTS – Rapid Transfer System

## Perfect results

- Specimen carrier loading in less than 5 secs
- Correlative LM/EM
- Rapid biopsy process

## Easy to install


- Factory mounted on EM PACT2
- Control via touch sensitive screen

## Easy to use

- Rapid loader for inserting specimen
- Automatic and reproducible specimen freezing

## Safe and convenient

- Outlet Dewar
- Safety lid
- Low noise
- Funnel for filling LN<sub>2</sub>
- Platform for cryopreparation, box, position for accessories



Schematic drawing showing the status of all parameters



Freezing data after ejection of the sample



# Leica EM RTS

## Application Solutions

**Specimen loading made simple with the Rapid Transfer System: insert specimen carrier into rapid loader. Automatically freeze by inserting rapid loader into RTS.**

Prepare your sample (eg leaf or a cell monolayer on a sapphire disc) ...

... to fit into one of the various specimen carriers (eg flat carrier  $\varnothing$  1.2 mm, depth 200  $\mu$ m or membrane carrier:  $\varnothing$  1.5 mm, depth 200  $\mu$ m). Specimen and carrier are held in the rapid loader. The long life pod is connected to the loading device with a bayonet lock.

Specimen and carrier are held in the rapid loader.

The long life pod is connected to the loading device with a bayonet lock.

The primed loading device then waits for freezing in the RTS.

By gently pushing the rapid loader into the RTS the carrier is tightened securely and then frozen automatically in less than 2.5 seconds.

After cryofixation the carriers are collected in the LN<sub>2</sub> bath of the EM PACT2 before freeze substitution in the Leica EM AFS2.





Leica Design by W. Hölbl

**The Tube Holder for every fluid you can think of ... blood, milk, cell suspensions ... nematodes – and more!**

Directly suck the suspension into copper tubes (inner  $\varnothing$  350  $\mu$ m) already mounted in their holder.

Alternatively, take up sample into cellulose microcapillary.

With the loading device the specimen can be inserted into the EM PACT.

After cryofixation the samples are stored under LN<sub>2</sub> in the transfer box.

Sample Punch: punching out of the central part of the specimen tube ready for Frozen Hydrated Sectioning or ...

... peeling away of the top of the copper tube for Freeze Substitution.

Storage of the copper tubes is carried out under LN<sub>2</sub>.

The specimen tube holders are recycled by reloading them with copper tubes.



# Leica EM RTS

## Application Solutions

### Taking microbiopsies for EM is now faster than ever before!

Prepare the Microbiopsy Transfer Station for RTS under the Optical Workstation...

... before taking a biopsy with the biopsy gun.

Insert the biopsy gun into the transfer station ...

... and transfer the tissue into the biopsy carrier.

The specimen and carrier are held in the rapid loader ...

... while the primed loading device waits for freezing in the RTS.

By gently pushing the rapid loader into the RTS the carrier is tightened securely and then frozen automatically in less than 2.5 seconds.

After cryofixation the carriers are collected in the LN<sub>2</sub> bath of the EM PACT2 before freeze substitution in the Leica EM AFS2 or cryosectioning in the EM FC6.







### Frozen Hydrated Sectioning made easy.

Alternatively, the carriers can be trimmed with a cutter ...



... for Frozen Hydrated Sectioning in the cryoultramicrotome.



### Preparation for follow-on procedures.

All necessary steps for follow-on procedures (eg. Freeze Substitution) can be conveniently performed in the cryopreparation box supplied with the Leica EM PACT2.



### The Freeze Fracture Holder for everything that can be fractured ...

Prepare the Freeze Fracture Station under the Optical Workstation ...



... and preload the freeze fracture carrier. A copper ring is put on top of the loaded carrier.



The carrier and copper ring are sandwiched securely in the pod with the supplied torque wrench.



With the loading device the specimen can be inserted into the instrument.



After cryofixation the carriers are collected in the LN<sub>2</sub> bath of the EM PACT2 before Freeze Fracturing.




Accessories are also available for loading EM PACT2 without EM RTS.




Mouse embryonic fibroblast  
grown on sapphire disc.  
Soazig Lelay, Jana Maentler and  
Paul Verkade, MPI-CBG, Dresden,  
Germany





Yeast frozen in the membrane carrier. Scale bar = 500 nm, 250 nm for insert.  
Courtesy of Mark van Breugel, Jana Maentler and Paul Verkade, MPI-CBG, Dresden, Germany.



Insulin granules producing INS-1 cells. Scale bar = 5  $\mu$ m.  
Courtesy of Joke Ouwendijk, Jana Maentler, Melanie Jäger, Michele Solimena and Paul Verkade, TUD and MPI-CBG, Dresden, Germany



Rat, Langerhans cells. Scale bar = 1  $\mu$ m  
Courtesy of Joke Ouwendijk, Jana Maentler, Melanie Jäger, Michele Solimena and Paul Verkade, TUD and MPI-CBG, Dresden, Germany

# Leica Microsystems – the brand for outstanding products

Leica Microsystems' mission is to be the world's first-choice provider of innovative solutions to our customers' needs for vision, measurement, lithography and analysis of microstructures.

Leica, the leading brand for microscopes and scientific instruments, developed from five brand names, all with a long tradition: Wild, Leitz, Reichert, Jung and Cambridge Instruments. Yet Leica symbolizes innovation as well as tradition.

## Leica Microsystems – an international company with a strong network of customer services

Australia:	Gladesville	Tel. +61 2 9879 9700	Fax +61 2 9817 8358
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Canada:	Richmond Hill/Ontario	Tel. +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Rueil-Malmaison	Tel. +33 1 473 285 85	Fax +33 1 473 285 86
Germany:	Bensheim	Tel. +49 6251 136 0	Fax +49 6251 136 155
Italy:	Milan	Tel. +39 0257 486.1	Fax +39 0257 40 3273
Japan:	Tokyo	Tel. +81 3 5435 9600	Fax +81 3 5435 9615
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Sollentuna	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Glattbrugg	Tel. +41 1 809 34 34	Fax +41 1 809 34 44
United Kingdom:	Milton Keynes	Tel. +44 1908 246 246	Fax +44 1908 609 992
USA:	Bannockburn/Illinois	Tel. +1 847 405 0123	Fax +1 847 405 0164

and representatives of Leica Microsystems  
in more than 100 countries.

The companies of the Leica Microsystems Group operate internationally in four business segments, where we rank with the market leaders.

### ● Microscopy Systems

Our expertise in microscopy is the basis for all our solutions for visualization, measurement and analysis of microstructures in life sciences and industry. With confocal laser technology and image analysis systems, we provide three-dimensional viewing facilities and offer new solutions for cytogenetics, pathology and materials sciences.

### ● Specimen Preparation


We provide comprehensive systems and services for clinical histo- and cytopathology applications, biomedical research and industrial quality assurance. Our product range includes instruments, systems and consumables for tissue infiltration and embedding, microtomes and cryostats as well as automated stainers and coverslippers.

### ● Medical Equipment

Innovative technologies in our surgical microscopes offer new therapeutic approaches in microsurgery.

### ● Semiconductor Equipment

Our automated, leading-edge measurement and inspection systems and our E-beam lithography systems make us the first choice supplier for semiconductor manufacturers all over the world.

 [www.em-preparation.com](http://www.em-preparation.com)

Leica Mikrosysteme GmbH  
Hernalser Hauptstrasse 219  
A-1170 Vienna Austria

Tel. +43 1 48899-235  
Fax +43 1 48899-350

  
**Leica**  
MICROSYSTEMS