




Leica M165 FC
Leica M205 FA / M205 FCA
User Manual

General Instructions

Safety concept

Before using your microscope for the first time, please read the "Safety concept" carefully. It is included in your instrument. It contains additional information about handling and care.



Bitte lesen Sie vor
der ersten
Benutzung
des Mikroskops
das Handbuchs
sorgfältig durch.
Es enthält
wichtige
Hinweise
zur
Bedienung
des Mikroskops
sowie
zur
Sicherheit.

Use in clean rooms

The Leica Micros can be used in clean rooms with humidity problems.

Cleaning

- Do not use any abrasive cleaning agents, chemicals or techniques for cleaning.
- Never use chemicals to clean coated surfaces or accessories with uncoated parts. This could damage the surfaces, and solvents could be carried into the coated parts.
- For most cases, we can provide special kits on request. Some products can be modified, and we can offer other accessories or use the existing ones.

Servicing

- Leica's may only be carried out by Leica Micros after the field service technicians. Only original Leica Micros parts may be used.

Responsibilities of person in charge of instrument

- Ensure that the Leica Micros are always operated, maintained and used by authorized and trained personnel only.



Important Safety Notes

User manual

The individual modules of the Jeeb M Series Technology Series include an interactive CD-ROM with relevant Jeeb Manuals in 20 other languages. Keep it in a safe place and read it carefully to the Jeeb M Series manuals and updates and also take care of you to download and install our website at www.jeeb.com.

This Jeeb Manual describes the basic functions of the individual modules of the Jeeb M Series Technology Series and contains important instructions for the operation of the Jeeb M Series and accessories.

The Safety Concept JeebSelfControl is additional safety information regarding the service work requirements and the handling of the Jeeb M Series Technology Series and electrical accessories as well as the safety instructions.

You can control the individual system services in the Jeeb SelfControl via the Jeeb M Series Technology Series. Please read the Jeeb Manual and the safety requirements in the Jeeb M Series.

Before installing, operating or using the Jeeb M Series, read the Jeeb Manual and ensure that you have the necessary safety instructions.

To install the Jeeb M Series, the conditions and the safety requirements of the Jeeb M Series Technology Series and the Jeeb M Series must be read and the safety instructions must be read.

Special manuals

Some manuals are provided for a number of accessories:

- Jeeb Manual M2-216-1 for the wiring of the high-pressure recovery system.
- Jeeb Manual for the 100-watt safety Jeeb M Series high-pressure recovery system.
- Jeeb Manual M2-267-1 for the installation of the system.
- Jeeb Manual M2-166-2 for Jeeb M Series.



Symbols Used

Warning of a danger

This symbol indicates a danger of electrical shock and/or fire.



The following can be harmful:

- Leads to damage
- Junction distances or degraded insulation

Warning of hazardous electrical voltage

This symbol indicates a danger of electrical shock and/or fire.



The following can be harmful:

- Leads to damage
- Junction distances or degraded insulation

Danger due to hot surface

This symbol indicates a danger of electrical shock and/or fire.



Important information

This symbol indicates additional information or instructions that need to be read carefully.



Explanatory notes

This symbol indicates additional information or instructions.

Figures

The numbers in parentheses within the descriptions relate to the figures and the items within those figures.




Safety Instructions

Description	Intended relations to the instrument or accessories and accessories to be used in the field of use	Responsibilities of person in charge of instrument
<p>The individual modules fulfil the highest requirements for observation and documentation of Leica stereo microscopes of their series.</p>	<p>Place of use</p> <ul style="list-style-type: none"> • Refer to 'Safety Concepts' booklet • Certain components must be placed at least 10 cm away from the eye and from flammable substances. • Avoid large temperature fluctuations, direct sunlight and vibrations. These conditions should start measurements and microscopic images. • The instrument and accessories are not suitable for use in the following circumstances: 	<ul style="list-style-type: none"> • Refer to 'Safety Concepts' booklet
<p>Intended use</p> <ul style="list-style-type: none"> • Refer to 'Safety Concepts' booklet 		<p>To be filled:</p> <ul style="list-style-type: none"> • The Modules, accessories and accessories are operated, maintained and repaired by authorized personnel in the field of use.
<p>Non-intended use</p> <ul style="list-style-type: none"> • Refer to 'Safety Concepts' booklet 		<ul style="list-style-type: none"> • Operators have read, understood and observed this Job Manual and carefully the safety instructions.
<p>The use of accessories, the components and accessories (e.g. on the eye) unless they are specifically intended for that purpose</p>		
<p>The instructions and accessories described in this operating manual have been tested for safety and reliable results. The responsible Leica affiliate must be consulted whenever the instrument is altered, modified or used in conjunction with non-ice components that are not intended for use of this manual.</p>		




Safety Instructions (Continued)

<ul style="list-style-type: none"> Repairs, service work Refer to 'Safety Concept' booklet Only original spare parts may be used. Before opening the instrument, switch off the power and plug in the power cable 	<ul style="list-style-type: none"> Integration in third-party products Refer to 'Safety Concept' booklet 	<p>This resort, especially the measures for reduction of the noise and vibration level set</p> <ul style="list-style-type: none"> On the arrangement of workspace, work assignments and work flow (changing tasks frequently). Through the training of the personnel, giving consideration to ergonomics and other safety aspects. The ergonomic design and construction of the case. If altered, it should be intended to reduce the level of the job-related strain.
<ul style="list-style-type: none"> Touching the live circuit can cause injury 	<ul style="list-style-type: none"> Disposal Refer to 'Safety Concept' booklet 	
<ul style="list-style-type: none"> Transport Before going to the assembly or repairing the individual modules of the case, if altered, it should be checked for necessary components In order to prevent damage, any violation should be avoided. If necessary, the case should be opened by the customer and checked then before use. 	<ul style="list-style-type: none"> Legal requirements Refer to 'Safety Concept' booklet 	
	<ul style="list-style-type: none"> EC Declaration of Conformity Refer to 'Safety Concept' booklet 	
	<p>Health risks</p> <p> Works with the altered IT case causes the following risks: <ul style="list-style-type: none"> Working in the case may lead to eye strain and injury to the eyes. Descending on the platform in the finished work situation and the work-related matters may occur. An </p>	




Safety Instructions (Continued)

 Direct contact with eyes can cause severe irritation or permanent eye damage.

The mask can be used to determine if any persons are wearing eye protection.

- Light sources: Safety regulations require measures of the structure.
- UV protection screen in front of the mask: The screen prevents the beam from being directly into the eyes.
- Dummy filter carriers in the free positions of the rapid filter changer prevent direct UV radiation from reaching the eyes.
- UV filters are installed in the observation beam path to protect the eyes.
- The safety grill protection on the mask housing prevents radiation from the mask.

Warning
 Radiation could damage the eyes. Therefore:

- Never look into the grill slot on the mask: This is a very bright UV radiation screen.
- Never look into the eyes: A reflection from the filter is in the beam path.
- Fill empty filter positions with dummy filter carriers (eg. V165 T1).
- Do not select a white, strongly reflective background in the mask: This



Safety Instructions (Continued)

Supply unit

- Always plug the supply unit from the same supply.
- When plugging and disconnecting the supply unit by
- Before plugging the supply unit by
- When replacing the high-pressure mercury lamp and other parts, such as the heat-absorbing filter or the collector
- During the maintenance work on the supply unit

Lamp housing

- Never open the lamp housing when the lamp is switched on. To do this, disconnect the lamp housing by
- Before opening the lamp housing, allow it to cool off for at least 15 minutes. Danger of explosion
- Never cover the conduct on the lamp housing. Danger of fire

Mercury lamp

- To do the user manual and safety instructions provided by the lamp manufacturer; and, in particular, the section on how to proceed in the event of a lamp break and release mercury.
- To prevent, remove the mercury lamp. To do this, it is advisable to use a vacuuming and sealed working device. The lamp housing is damaged by the transition elements.
- When the required level of isolation is followed, the manufacturer's specifications and the rules refer to the necessary unit
- To do this, the necessary isolation, as well as decontaminated mercury, is necessary
- These measures ensure that there is no damage caused by emitting, non-toxicly isolated or reusable used mercury lamps



Contents

General Instructions	2	leica_m_series - Preset Lists	30
Mounting Society Vials	3	leica_m_series - Assembly	31
Synopsis Used	7	leica_m_series - Adjusting the Resolution	33
Society Vials Lists	5	leica_m_series - Centering the Objectives	35
		leica_m_series - Terminal Jit Protection Screen	36
The Leica M Series for Fluorescence		Terminal Head - gnt Base set up - 370; Before Terminal Job	37
Configuration	13	Terminal Head - gnt Base set up - 370*	38
A Step Towards Infinity	17	set up - VC; set up - VC	39
The Electronics Control, Cornerstone and Society of Your Case - Items 15	15	set up 3070 Menu - Mechanics Stage Assembly	40
The Module - Design: Everything is Relative	16	set up 3070 Motor and Mechanics Stage Assembly	43
Motor Control - by	17	Series Connections	47
On MeCa	18	Series Case Data	48
		Series Degree	49
Assembly		set up 225000 VC	50
Assembly by the "ocusing Camera" - Jobs	20	set up 225000 VC: Alternative Assembly	51
Assembly by the "ocusing Camera" When "xident" gnt Base	21	set up 225000 T: Assembly	52
Motor and "ocus Test" of the "ove Belt"	22	set up 26000 - Society Vials	57
Optics Center	23	set up 26000 - Preset Lists	57
Table	27	set up 26000 - Connection to the "Junction" Assembly	60
Types etc.	25	set up 26000 - Test of the Series	62
Objective	26		
Jit Protection Screen Assembly	27		
Objective Vials etc - Assembly	28		
Objective Vials etc - Adjusting Resolution	29		

Contents (Continued)

Quick Start Guide	
The Fastest Route to Success	67
Overview of the M Series Microscope	65
The Correct Filter Setup for D-SLM	66
Buying the Dyes	67
Working	68
Adjusting the Test Slide on the T.O.D Drive	69
Changing Magnification (Zoom)	70
Key Steps and Magnification Levels	71
Recovery: More Control and Convenience on Your Work- station	72
	73
Eye Pieces	
Magnification Factors of the Eyepieces	75
Interfiliers	76
Correct Correction	77
Correct Correction and Reliance	78
Calculations	79
Photography & Video	
Photography & Video	81
Photography and Components	82
Photography and Video 50%	83
Photography and Video 100%	87
Photography and Video 100%: Digital Image 3-75	85
Microscope Carrier	
The Objective Noses	87
Objectives and Optical Accessories	
The Different Types of Objectives	89
Slides	
Slide - S Series Head - gnl Slide Carriage	91
Slide - S Series Head - gnl Slide Carriage	92
Slide - S Series Head - gnl Slide Carriage	93
Slide - S Series Head - gnl Slide Carriage	94
Slide - S Series Head - gnl Slide Carriage	95
Slide - S Series Head - gnl Slide Carriage	96
Leica TL RCL The Deflection Mirror	97
Slide - TL: Carriage Interfiliers and Temperature	98
Slide - TL: Slide - TL: Objectives	99
Slide - TL: Methods of Test Head - gnl	100
Slide - TL: Test Results	101
Slide - TL: Test Results	102
Slide - TL: Test Results	103
Slide - TL: Test Results	104
Slide - TL: Test Results	105

Contents (Continued)

System Illumination	
Leica C35000 MC	107
Leica C35000 TL	109
Leica L6000 – About the Instrument	110
Leica L6000 – Operation	112
Leica L6000 – Troubleshooting	113
Centering and Meridian	117
Fluorescence	
Filter Change	116
Filter Fluorescence Filterly Manager	117
About Fluorescence Microscopy (Dry Assay) for Leica M205 TM	118
Load Filter Change and Filter Types	119
Slide Filter Holders	120
Switching the Filter Change	121
Observation without Fluorescence	122
Double Exposure	123
Comments on the Fluorescence System	127
Accessories	
Leica Hand Wheel and Tool Set Kit	126
Leica Slide Move	127
Leica Slide Catch	128
Leica PSC Controller	129
Dimensional Drawings	
Leica M165 TC	131
Leica M205 TM	132
Leica M205 TCA	137
Leica C35000 MC	140
Leica C35000 TL	141
Leica M165 TC	142
Specifications	
Leica M165 TC / Leica M205 TM / Leica M205 TCA	142
Leica EL6000 – Specifications	143
Leica L5000 – Standard Method	150
Leica L3000 – Standard Method	151
Leica L1000 – TC	152
Leica L5000 Diga	153
Leica 3000 Motorized XY Stage	157
Appendix	
Calculating the Total Magnification and Field of View Diameter	156
Cell, Microtome, Carrier, Coverslips	157

The Leica M Series for Fluorescence



Congratulations!

Congratulations on completing your new Leica M Series selected microscope. We've worked to exceed your expectations as never before, having learned our decades of experience in the areas of optical technology and ergonomics to create a microscope that works by the way.

The Leica M Series addresses the requirements of the new microscope systems: excellent objectives, high-quality eyepieces, and, by the way, the modular design ensures that the M Series adapts perfectly to your needs – via the optional accessories you require in your lab.

The new imaging system, including the smart, objective and eyepiece, is apochromatically corrected with much technological effort. Contrast, sharpness, richness in detail, resolution, image and color fidelity are optimized. In addition, the selected illumination system delivers at every zoom level that light utilization is at a maximum and that fluorescence images are intensely illuminated, set against background.

Though the Leica M Series addresses all selected microscope systems, especially, when you select the Leica M Series, you require a certain degree of care and attention. Therefore, we recommend that you read this manual. Learn to use the M Series as you need it, regard it as a tool, and the M Series will serve you best. The M Series will be your selected microscope – can't you to work as smoothly and reliably as on the very first day.

We wish you the best of success in your work – please, you are now equipped with the best tool.



A Step Towards Infinity

Ever since the introduction by Harold S. Greenough, stereo microscopes have worked according to the optical principles devised primarily in Ernst Abbe's research. For over a century, ingenious optical designers and engineers have worked to push magnification, resolution and image fidelity to the limit permitted by optics.

Today so, they have kept seeing what is hidden by the interaction between three factors: the angle a microscope's resolution, the size the objective works in, and the waves the diameter of the objective. The thread remains the same: the objective becomes distorted – a sphere becomes an ellipse, a flat surface curves towards the objective.

How do you get into the microscope and how do you get out?

Limits are made to be broken.

The Leica M205 FA is the world's first fluorescence stereo microscope with a four-engine design of 20.5X. This accounts for the increase in resolution of the microscope. With the new design of the Leica M205 FA, the microscope has exceeded the resolution of the microscope beyond the 1.050 μm limit. The increase in magnification, the resolution, which corresponds to a resolution of 0.72 μm .

Of course, this increase in performance benefits your everyday work. Set up your microscope on the microscope like with computer-aided design of microscope and discover details in stereo microscopy that you could never see before.




The Electronics: Comfort, Convenience and Safety for Your Experiments

Never before have electronics been used as extensively in science as in the past. With series, parallel, series-parallel, and other circuits connected to your electronic controls – which provide a number of advantages.



Control is not only transmitted, but also supply the power.

Reliability for your experiments. Select a wide variety of results and increasing your command, accuracy, precision and development. The continuous recording captures parameters such as the magnitude of the current, the position of the



LISA Application Suite (LISA) evaluates the data collected to and can solve the test situation.

designer and more end users in the information in the case. AS software thus you always know the conditions under which an image was acquired.

Power cables

A significant part of the cables have been routed in the interior of the cabinet. Delete also the cables via the interface between the control and the control center and the control. The direct result is that you need fewer cables – this not only saves your workstation space and more comfortable, even makes it easier to get.



The Modular Design: Everything is Relative

The Leica M series provides maximum flexibility through equipment, lenses and the modular configuration and the compatibility that Leica has established for decades. The optical center, eyes, ears and nose can be configured in any way you choose, allowing you to create the microscope that best suits your needs.

Just like this, you will notice that the components and individual components do not differ significantly. Whichever configuration you choose, you will quickly see digital nature.

Have a special request? Let us know! We will explore every possibility and let you know when it comes to developing custom-specific solutions. If you have a special request that will be met, please visit us at www.leica.com. We have a solution for every project.



Maximum Compatibility

Jetco eyepieces were created to ensure that the new Jetco M-series microscope objectives could be carried over in existing series. This means that objective lenses, tubes and so on can be reused.

Objectives

All new objectives of the Jetco M-series are available. Nothing that when used with the objective tubes etc. They can be used with existing M-series then tube objectives.

If you prefer, you can still use the previous Jetco objective series. In this case, the objective is no longer available.

Tubes

The difference between the old series and the new series is the new series of existing tubes fit the new M-series. The new tubes are designed for eyepieces with field number 23, while the predecessor models were only designed for field number 21, resulting in a smaller object field.

Eyepieces

However, the new M-series eyepieces have an eye ring and long clicks to provide immediate feedback of accurate adjustment.



On We Go


If your first jetted exercise has already been prescribed and confirmed by your jetted consultant, click on the table below to learn more about jetted exercise and go directly to the Quick Start guide on [page 63](#).




If, on the other hand, you are entering the "exercise your jetted consultant" category, click on the table below to learn more about jetted exercise and go directly to the Quick Start guide on [page 19](#).



If everything you need to start about the correct use of fluorescence-related parts, refer to [page 115](#).



Assembly



Assembling the Focusing Column for TL Bases

The first step is to connect the focusing column

of the M series to the corresponding base

Assembling the column adapter

1. Securely fasten the column adapter on the column by the four torx screws



Assembling the focusing column

2. Securely screw the focusing column to the base by the four torx screws



Tools used

- Torx set screws: set, 3 mm



Assembling the Focusing Column With an Incident-light Base

When using an incident-light base, the focusing column and motorized focus are tableted directly on the base; no external tableting is required.

Tools used

- Torx screwdriver, 3 mm

Assembly

1. Place the focusing column on the tableting




3. Screw the base securely onto the focusing drive




2. Insert the four screws provided into the guide holes of the base



Motorized Focus: Restricting the Travel Path

 Depending on the work situation, it is best to restrict the travel path of the motorized focus. This prevents injury during handling of the microscope. By fingers or hands becoming trapped or the specimen being damaged, even being damaged by it.

Readjusting the motorized focus
The motorized focus is factory adjusted and the yokes will need to be readjusted—even if the travel path has not changed.

 **Caution:** The barrier is above the motorized focus. Moving the position deck or dist. This case, the specimen that is released using the motorized focus near the Specimen Touch Pad will cause contact the specimen table.

Restricting the bottom travel range

1. Move the motorized focus into the lowest position you will use.
2. Loosen the screw of the travel stop on the side of the focus yoke arm.
3. Push the travel stop to the height of the motorized focus.
4. Loosen the screw of the travel stop by moving the screwdriver handle and moving it upwards.
5. Tighten the screw of the travel stop.



Optics Carrier

Tools used

- Torx screwdriver, T10

Assembling the optics carrier

1. Place the optics carrier on the focusing column so that the screw fits into the thread provided and the lug fits into the groove.



2. Press the optics carrier downwards to the focusing column screw. Tightening your other hand.



Tube

All intermediate tubes that fit between the optics carrier and the binocular tube are fitted in the same manner:

Tools used

- No tools required

Preparations

1. Fix the positioning screw and remove the protective cover:



Assembling the tube

2. Push the tube for exercise, the fixed screw is pushed into the dovetailing end while its other end fits into the guide groove of the positioning screw. Press it with the guide groove.



3. When sliding the tube only slightly, gently tighten the positioning screw. This adjustment is brought to the correct position.



Eyepieces

Tools used

- Your hands required

Magnification range

You can extend the overall magnification range by using 10X, 16X, 25X and 30X wide-field eyepieces for persons wearing glasses.

Preparation

1. If you will use an ocular guard, refer to [page 79](#).
2. Remove the ocular guard.



Inserting the eyepieces

3. Push the eyepieces into the tubes as far as they will go and check to ensure that they fit tightly and accurately.



4. Securely tighten the locking screws.



Objective

Tools used

- No tools required

Preparation

1. Remove the objective caps on the objective lenses by turning 1.



Attaching the objective

- Hold the objective firmly during assembly and disassembly so that it does not slide onto the stage plate. It is also a delicate piece for the 2X objective, which is delicate, which is very heavy. Therefore, select the objective stage plate first.

2. Screw the objective clockwise into the objective turret.



Alternative fastening options

- Using the objective bases see [page 28](#) instructions on [page 28](#).
- Using the objective adapters see [page 30](#) instructions on [page 30](#).



UV Protection Screen: Assembly

Tools used

- Allen key


Use with care objective


1. Adjust the UV protection screen elements using the screw.
2. Tighten the hexagon-head screw.
3. Adjust the UV protection screen using the screw.
4. Tighten the hexagon-head screw.
5. Tighten the UV protection screen using a hexagon-head screw to the left-hand side of the "cross-axle" etc.

Intended use

The UV protection screen in front of the backrest prevents the operator from looking directly into the UV rays.

Safety Notes

 UV radiation can damage the eyes. Therefore, it is mandatory to make the UV filter and adjust it correctly.

 Always position the UV protection screen so that the operator can never look directly at the grill seat.



Objective Nosepiece – Assembly

Preparations

Hold the objectives firmly during assembly and assembly so that they do not fall out of the stage plate.

2. Turn the rotating set by 90° and align the objective tubes to the drive housing for the front. 1. Screw the objective nosepiece firmly into place.



- Move the drive housing, the key strands and remove the oil cover. The center hole is ready for use.

Assembly

1. Turn the front part clockwise for the objective tubes etc.



3. Insert the three 21 screws on the objective front of the oil cover and remove the intermediate ring.



4. Screw the oil cover onto the objective tubes etc.

5. Screw both objectives onto the objective nosepiece. It makes no difference which side of the objective tubes occurs.

6. Insert the locking screws on both sides of the objective tubes etc.



You can't work at the surface by the rotations on the next page.



Objective Nosepiece – Adjusting Parfocality

The following procedure may be necessary to achieve parfocality between both objective lenses in a microscope. Parfocality means that the specimen remains in focus when the objective lenses are swapped.

This procedure may be needed if you receive either of the two objective lenses in the table:

The following steps ensure the parfocality of the 1X and 2X objective lenses. You are using either objective camera field of view or the 2X objective in the description with the objective with the stronger magnification.

Preparation

- Use the standard eyepiece.
- Set the dust cap correctly on the eyepieces.


Adjustment

1. Rotate the 2X objective into the center field and set it to the lowest magnification.
2. Focus on the specimen.
3. Rotate the 1X objective into the center field.
4. Turn the objective on the thread in both directions until the specimen appears sharp.
5. Toggle to the 2X objective.
6. Select the strongest magnification and refocus until the specimen appears sharp.

7. Toggle to the 1X objective.
8. Turn the objective on the thread in both directions until the specimen appears sharp again.
9. Repeat steps 4 through 6 for the series of objective lenses. Repeat the checks with the other objective lenses as well before you set the procedure.
9. Tighten the locking screws.



Leica FluoCombi III – Preparations

 Before setting up the microscope, you have to disconnect the instruments for the camera and the microscope. These instructions can be found in the manual and the connected instruments.

Preparations

1. Turn the two screws used as the instrument holders.



2. Screw the two screws into the holes provided at the stage.



In this way, you can use the microscope as the instrument holder if necessary.

3. Screw the SX objective into the microscope.

4. Screw the other objective into the objective holder.



Leica FluoCombi III – Assembly

1. Guide the **JoCombs** into the **seal**. While holding it flat, ensure that it is aligned with the **marker** on the **focus ring** at the **end** of the **grill**.
2. Hold the **JoCombs** secure and screw in the **focus ring** as far as the **provided screw**.



Leica FluoCombi III – Assembly (Continued)

Assembling the optics carrier

1. Fix the microscope and the objective on the optics carrier:



2. Align the optics carrier to the microscope body and screw it into place.

3. Open the side integration cable.




4. Assemble the tube (see page 24).
5. Assemble the objective (see page 50 ff.)



Leica FluoCombi III – Adjusting the Parfocality

Setting up the objective nosepiece

 The following configuration only has to be carried out once. It ensures that the objective lenses are aligned and that the distance does not need to be brought into focus every time.

1. Place base then under the microscope
2. Take the selected electronic objective under the base unit
3. Adjust the view so that you can see the side then the over-view.
4. Center the base then under the base unit
5. Adjust the magnification to the highest level and focus on a flat specimen.
6. Take the 5X selected under the base unit

7. Adjust the carrying screw on the objective



8. Adjust the view so that you can see the side then the over-view.



Leica FluoCombi III – Adjusting the Parfocality (Continued)

9. By turning the objective ring, focus on the slide. Then tilt the slide. The slide should be adjustable in 5°.





11. Tilt the slide.
12. Tighten the centering screws. (See Step 7).
13. Tilt the slide to a different objective under the same magnification. This causes the two objectives to be out of focus. This means that when the magnification is between the two objectives, the slide is out of focus.

10. Change the magnification for the center of the specimen to be visible in the field of view.



Leica FluoCombi III – Centering the Objectives


 The following configuration only has to be carried out once. It ensures that the specimen remains in the optical center when the objective is changed.

 To allow fast switching results in the following settings, insert the graticule into the eyepiece. This ensures that the center stays constant.

Centering the objective

1. Take the 5X objective under the specimen.
2. Focus on the top of the specimen.
3. Take the Leica 5X electronic objective under the specimen.

4. Center the objective over the specimen by observing and lighting the specimen on one side of the objective handle while turning the specimen handle in the opposite direction.



Leica FluoCombi III – Filter and UV Protection Screen

Inserting the dichroic mirror


1. Pull the dummy filter carrier out of the holder:



2. Insert the correct sized rigid carrier into:

Assembling the UV protection screen

1. Screw the UV protection screen to the holder provided on the objective carrier:




2. Orient the UV protection screen so the you cannot see directly into the objective.



Transmitted-light Base Leica TL BFDf: Before First Use

Removing the transport anchors

⚠ Before you can use the transmitted-light base for the first time, it is absolutely necessary to remove the two transport anchors.



Transmitted-light Base Leica TL BFDf

Standard delivery

The base is delivered with the objective slide installed. The selected slide (see 3020a microscope slide or Leica slide) slide (10450562), and the focusing drive will have to be mounted etc.

Stage assembly

The Leica T-370 Transillight will be delivered equipped with three different stages.

- Leica 3020a Transillight microscope slide
- Leica 3020a Transillight microscope slide
- Leica slide disk 10450562

The selected stage is mounted on the base before continuing. You can switch between the stages of any Transillight microscope slide.

The focusing stage is assembled on the base with the stage mounted. To assemble is performed in reverse order of the drawing.

Standard stage

1. Turn the glass plate over the microscope glass microscope slide.
2. Push on the stage on the Transillight with the base. Push the four screws on over those of the base.
3. Push on the stage to the base with the four screws. The stage need screws.
4. Turn the glass plate back into the microscope slide.



Leica TL RC / Leica TL RCI

The base is delivered with the edelcr: side visible. The selected stage is also shown in the reference stage on the side of stage 10450562), and the focusing drive will have to be mounted etc.

Ensure that the instruments are unobscured and view, adequately sized and mounted under eye.

Stage assembly

The stage TL RC / TL RCI is a 100 mm diameter stage equipped with three different Leica stages. The selected stage is mounted on the base before continuing with the assembly between the stages of the TL RC / TL RCI and the focusing drive.

The focusing stage is not assembled on the base without the stage mounted. Otherwise, it will be damaged. The order of the focusing stage is as follows:

Standard stage

1. Take the glass plate from the reference stage on the side of stage 10450562).
2. Push the stage on the TL RC / TL RCI over the base of the base.
3. Align the stage to the base on the TL RC / TL RCI.
4. Push the glass plate back into the side of stage.



Leica IsoPro Manual Mechanical Stage: Assembly

Leica IsoPro Mechanical Stage


Before the steel sledge mechanism is mounted to the base, the control buttons are attached to the left and right sides of the mechanism sledge.

The controls are to be mounted on the left-hand side, the gear rod on the right-hand side of the mechanism sledge must be inserted and reflected in reverse.

1. Take the gear rod from the mechanism sledge.
2. Turn the mechanism sledge around and place the gear rod's surface.
3. Change the gear rod from the left to the right-hand side.
4. Set the next two steps to take the controls.

Before the steel sledge mechanism is mounted to the base, the control buttons are attached to the left and right sides of the mechanism sledge.

Left-hand operation
The controls are to be mounted on the left-hand side, the gear rod on the right-hand side of the mechanism sledge must be inserted and reflected in reverse.



Leica IsoPro Manual Mechanical Stage: Assembly (Continued)

- | | Control assembly | Mechanical stage assembly |
|--|---|---|
| 1. Take the glass plate from the reference stage and turn it around. | 1. Take the glass plate from the reference stage and turn it around. | 1. Place the reference stage on the base. |
| 2. Change the gear and from the left to the right-hand side. | 2. Align the axis in the control stations to the desired reference axis of the reference stage (reference X). | 2. The user sets the reference stage carefully towards the user, using the lower side of the reference stage. |
| | 3. Align the axis in the measurement head (reference Y). | |
| | 4. Align the cover plate to the reference stage. | |



Leica IsoPro Manual Mechanical Stage: Assembly (Continued)

3. Affix the mechanical stage evenly to the three needed rods.
 1. Insert the extension plate into the base using the Allen key provided.
 2. Affix your focusing drive column using screws using the three hexagon-head screws.
 3. Reflect the extender plate to its original position using the three hexagon-head screws.
4. Now, move the mechanical stages' axes to your original location of the column.
5. Insert the glass plate into the mechanical stage.




Leica IsoPro Motorized Mechanical Stage: Assembly

Basics

The motorized grinders of the Leica IL series (IL 8HD, IL RC, IL RC0) are supplied with a table ed extension. These extend stage Leica IsoPro mechanical stage or standard stage (0.450.562) and the focusing drive will have to be mounted etc.


The motorized mechanical stage is suitable for precision adjustment. During table start avoid subjecting the stage to mechanical severe vibrations.

1. Check the mechanical stage for the standard accessory and install on the motorized grinder.



1. Microscope base
2. Motorized mechanical stage

2. Secure the motorized mechanical stage to the base using three W screws.






3. W screws



Leica IsoPro Motorized Mechanical Stage: Assembly (Continued)

3. Remove the two M3 screws and the handle from the mechanical stage.
4. Remove the two W screws and the sleeve from the mechanical stage.
5. Remove the four stainless-steel caps from the mechanical stage.



 After removal, keep the stainless-steel caps in the plastic bag provided for future reference.




Leica IsoPro Motorized Mechanical Stage: Assembly (Continued)

The bases

The IsoPro motorized mechanical bases can be equipped with three different stages: Standard stage, the standard stage, the IsoPro mechanical stage. The selected stage is mounted on the base before carrying out the assembly. You can select between the stages at any time with the external movement.

The driving direction assumed for the base without the stage mounted. Otherwise, a warning message order of the driving

Mechanical stage and base

 Never move the sledge of the motorized mechanical stage in the X direction, as otherwise the mechanical system will be damaged.

1. Place the mechanical stage onto the base.
2. Push the upper part of the mechanical stage carefully towards the bearing end, when the cover is onto the IsoPro filled-grip base.
3. Attach the mechanical stage evenly to the three threaded rods.
4. Now, move the mechanical sledge as far as you go in the direction of the column.
5. Attach the glass plate to the mechanical sledge.



Leica IsoPro Motorized Mechanical Stage: Assembly (Continued)

Mechanical stage for X-Y stage control module

1. Plug the 2-pin X-Y stage control module into one of the free connector ports.

2. Plug the 2-pin X-Y stage control module into the 2-pin interface

- The 15-pin D-sub interface is also available for connecting other controllers, such as the Leica foot switch (10 447 398) or IL KLI 1200 (10 447 352).

- The 15-pin D-sub interface is provided for use with the Leica S-Move controller (11 501 197).

3. Plug the sensor cable into the socket provided and into a grounded sensor socket.

As soon as the adjustments have been jogged into the control box and they have been jogged into the sensor supply, the mechanical stage will be automatically driven to the end position.

4. Connect the X-Y stage cable and X-Y Stage I/O Module (using a suitable USB cable).


Cables: Connections

The new case M series features extensive encoding which allows for accurate date and serial tags can be read out from either the AC end or audio end.

- The interface to the audio connector is on the bottom



The terminals



The connection to the AC end to other M series terminals is made using the terminals on the rear side of the connector.

1. 3-pin terminal for auxiliary equipment such as the ring amplifier; the base end other accessories with the case provided cable
2. Terminal for the power supply provided
3. USB terminal for the connection to the AC.



Cables: Cable Duct

The integrated cable duct in the carrying case is a mesh case that surrounds the microphone. For example, the JS3 or TWS versions of the case are not designed to use a duct.

Feeding the cables

1. Loosen the three screws on the cable duct.



2. Remove the cover of the cable duct.

3. Place the cables in the cable duct and secure the cover as originally.



3. Do not tie the length of the cable ends you need before screwing on the cover. For thick cables, it is difficult to change the length retroactively.



Leica LED5000 MCI

The Leica LED5000 MCI for Multi Contrast JET is not suitable as light source for all microscope objectives. The objective should be removed during assembly, see page 231.




Constraints

The Leica LED5000 MCI cannot be used together with the objective 'Hofmann' etc.

Assembly

1. Hold the LED5000 MCI in one hand and lightly screw the objective into its position as described in the manual.



2. Connect the cables as provided to either of the two sockets. (The flat part of the plug must be facing downwards.)



3. Plug the other end of the cable 'Hofmann' into the 'Hofmann' sockets on the LED JET.



Leica LED5000 MCI: Alternative Assembly

Under certain circumstances, the grill source that will be moved along with the optical center of the microscope is not a good option, in which the Z-axis changes with the angle of rotation of the grill table. In the case of such a situation, the Leica LED5000 MCI is directly installed to the column.

Installation on the column

1. Pull the telescope out of the Leica LED5000 MCI.



2. Screw the telescope into the column at the top using the following screws. You tighten the screws; the telescope slides a certain distance toward the correct position.




3. Push the Leica LED5000 MCI onto the telescope.



Leica LED5000 RL: Assembly

Required tools

- None

 The Leica LED5000 RL is installed on the objective using a 3-way screw. This screw can be used to secure the distance between 30 mm and 80 mm.

Constraints

The Leica LED5000 RL can be used only in conjunction with the telecentric 1X and telecentric 0.62X objectives. With the other objectives, the working distance is too short or adequate for telecentric.

The 3-way adapter cannot be used together with the objective microscope.

Assembly

1. Connect the Leica LED5000 RL as described in the 3-way adapter.



2. Place the Leica LED5000 RL against the objective from below, push it into place as shown and screw it into place.




3. Plug the other end of the cable into the external USB socket on the LED5000 RL.




Continued on next page.







Leica LED5000 RL: Assembly (Continued)

 If you work with the remote focusing camera without integrated electronics, the Leica LED5000 RL must be supplied with some external power supply unit that is included in the delivery package. Moreover, the controller must be controlled via the application software in this case.

 To get full accessibility of the software, the original controller should be installed with the cable using accessories. It is also possible to turn the original controller's display or external display unit on by using the Leica LED5000 MC system controller. Moreover, the original controller must be connected to the Leica LED5000 MC.





Leica EL6000 – Safety Notes


-  The manufacturer assumes no responsibility for any use outside of the intended use or use outside of the specifications for Leica Microsystems Wetzlar GmbH or any resulting instruments.
- In all cases, the Declaration of Conformity shall be valid.
-  This IVD device is not intended for use in the patient environment defined by IEC 60601-1-1. It is also not intended to be combined with medical devices as defined by IEC 60601-1. If a microscope is electrically connected to a Leica instrument, accessories with IEC 60601-1, the requirements defined in IEC 60601-1-1 shall apply.
- General safety notes**
 This Leica EL6000 device has not been tested in accordance with:
 EN 61010-2-101:2002
 EN 61010-1:2001,
 EN 1010-1:2001.
 Safety requirements for electric equipment for measurement control and regulatory use.
- The devices and accessories described in this operating manual have been tested for safety and suitable results.
- The responsible Leica affiliate or the main plant of Leica must be consulted whenever the device is altered, modified or used in conjunction with other components that are not described in this manual.
-  The manufacturer is liable for the safety of the instrument. All accessories will be tested and certified in accordance with:
 EN 61010-2-101:2002
 EN 61010-1:2001,
 EN 1010-1:2001.
-  Safety safety defects or malfunctions of the Leica EL6000 are identified, the instrument must be returned to the manufacturer for the service system and secured against additional use. The Leica EL6000 must be sent to the manufacturer or authorized representative of the manufacturer.





Leica EL6000 – Safety Notes (Continued)


 The Leica EL6000 connected grid source generates high energy grid voltages. Accidents there can also be triggered or caused by grid level out directly into the grid guide output of the instrument or into the output of the grid guide connected to the instrument. Always ensure that the output of the connected grid guide is securely connected to the system to be grid level before starting on the connected grid source. The maximum voltage is 100–270 V AC, 50–60 Hz. With this voltage and frequency range, the instrument always reacts to the connected power supply. Observe the instrument's maximum supply voltage and do not destroy the instrument and the connected components.

 To prevent damage to the grid guide connected to the device, the heat-absorbing filter inserted into the device must be free of damage (e.g., cracks etc.). The heat-absorbing filter is not inserted or is defective, the service life of the grid guide will be decreased.

 The use of the error meter in the instrument can lead to safety hazards. At the end of the service life, based on the error according to the service regulations, the meter's accuracy of the error provides corresponding information, or otherwise the error is stated in the documents. To avoid the manufacturer's liability, please use the error meter.

 The Leica EL6000 connected grid source is designed for a voltage range of 100–270 V AC, 50–60 Hz. With this voltage and frequency range, the instrument always reacts to the connected power supply. Observe the instrument's maximum supply voltage and do not destroy the instrument and the connected components.

 Only uses of the specified type and the stated current may be used as accessories. The use of incorrect tools may result in an accident.

 The instrument is intended exclusively for operation in dry areas. Do not use the instrument outdoors or in areas with high humidity.



Leica EL6000 – Safety Notes (Continued)

Electrical safety	
For wear device	
Supply voltage	100–240 V - ±10 %
Frequency	50–60 Hz
Power consumption	Max. 210 W
Inputs	5x20, 2.5 A, 3-wire wiring system, category II, see page 117
Ambient temperature	0–40 °C
Relative humidity	10 %–90 %
Overvoltage category	
Protection degree	2



Leica EL6000 – Preparations

Setting up the instrument

⚠ Set up the Leica EL6000 so that the front is easily accessible and visible.

- The ventilation slots on the sides and back side of the device must not be covered.

- Make the device warm up at least 150 minutes before use of the instrument.

Instructions for lamp replacement

⚠ The EL6000 carries a mercury lamp. Be especially careful to follow the handling instructions and safety rules provided with the EL6000.

⚠ Make the device inaccessible when using the EL6000. Remove the EL6000 according to the instructions.

⚠ The EL6000 is used in the instrument and becomes very hot during operation and may reach a high temperature when not before changing the EL6000. It is mandatory to leave the instrument to cool off for at least 20 minutes.

⚠ The Leica EL6000 compact light source comes with a factory-adjusted EL6000. The EL6000 should not be used without adjustment. This may cause the danger of damage to the EL6000 during transport.

⚠ Before transporting the Leica EL6000 compact light source, you must have to remove the EL6000.

Inserting the lamp

1. Put the cover bag out of the socket so that the EL6000 is disconnected from the power system.
2. Fix the cover with screws of the cover using a 3 mm Allen key.



3. Remove the top cover.

Leica EL6000 – Preparations (Continued)

4. Carefully pry the Leica EL6000 on its side and then open by loosening the screws.



5. This view is used to label the air.

5. Push the pressure seal back towards the front to clear the eye.



6. Label the air.

7. To label the groove in the connector, label the air. Mark the air in the corresponding tag in the air tag.



Leica EL6000 – Preparations (Continued)

7. When the ears are fully connected in the ears mount, release the ears:

8. Connect the plug of the ears with the cords of the instrument.



9. Check to see whether the ear-mounting filter has been inserted into the shaft in order to protect the connected grid guide.




Ensure that no wires are touching the reflector of the ears.

10. Close the housing cover and tighten the screws.



Leica EL6000 – Connection to the Fluorescence Lamp Housing

 Always connect the grill guide to the emission housing of the selected microscope to prevent damage to the user with the high-energy grill.

1. Insert the isolation screw and remove the protective cover of the fluorescence emission housing.



2. Insert the adapter into the fluorescence emission housing and tighten the locking screw.



3. Insert the threaded protective cap on the grill guide.




4. Insert the start end of the grill guide into the emission housing of the selected microscope as described in the guide.

5. Securely tighten the emitting screws.




Leica EL6000 – Connection to the Fluorescence Lamp Housing (Continued)

6. Insert the original or the grid guide into the grid outlet of the Leica EL6000. There must be no play at all.

 Only use grid guides: the grid must be of the type 'Starz ang', as otherwise damage to the instrument and danger to the user can result; see also page 22.



 Before you open the shutter, the grid guide must be connected on both sides. Otherwise, the rotating grid can cause injury to eyes and skin and damage to the 'Live View' box. Insert the grid correctly into the grid guide.

7. Connect the Leica EL6000 to the zoom system by the zoom cable.



Leica EL6000 – Replacing the Lamp

The procedure for replacing the lamp does not differ from the installation procedure outlined on [page 57](#), except for the fact that the starting lamp has to be removed first.

Replacing the Lamp

The air is cooled in the instrument becomes very hot during operation and the high static pressure when hot. Before changing the air, the mandatory task on the instrument to cool off for at least 20 minutes.

1. Put the cover bag out of the vessel, so that the EL6000 is disconnected from the power system.
2. Fix the 20 mm screws of the cover using a 3 mm allen key.

3. Remove the housing cover.
4. Carefully by the Leica EL6000 on both sides that the heating coils stay.

This makes it easier to take the air.

5. Put the pressure ball back towards the front side using the lever.
6. Pull out the air and disconnect the bag from the cooling in the instrument.
7. Take the new air as outlined on [page 58](#), by using the lever.



Quick Start Guide




The Fastest Route to Success

Your Leica Stereo microscope has been delivered, unpacked, ready for use. To get the best condition for your Leica Stereo microscope, you will need to get it into working condition. Your next step should be to study the Quick Start Guide without missing the most important steps and details.

This manual will then familiarize you with the finer details of your microscope. The following pages contain important, precise information that needs to be followed every day:


• Getting the best from your instrument



The Microscope Quick Start Guide (included with the instrument)

Overview of an M Series Microscope

1. Eyepieces for spectacle wearers with diopter correction and eyecups
2. Trinocular tube
3. Fluorescence housing
4. Filter Changer
5. Display
6. Motorized Focus
7. Objective No. display
8. Lens Shift Touch
9. MP photo start screen



The Correct Interpupillary Distance

The interpupillary distance is correct if you see a single circular image field when looking through them.

If you see a double image field, you may need to adjust the interpupillary distance. Adjusting the interpupillary distance will allow you to see a single image field.

Reference value

The distance between eye and eyepiece measures approx. 22 mm for 10/238 wide-field eyes and approx. 19 g. eyes.

Adjusting the interpupillary distance

1. Look into the eyes.
2. Hold the eyes at the same height. Just the eyes, not the eyepieces. If you see a double image field, you need to adjust the interpupillary distance.
3. Slowly approach the eyes. As you approach, you will see the image field without coming cutting.



Using the Eyepieces

The eyepieces form the connection between the tube and the eye of the observer. They may also form a part of the tube and they are ready to use.

Each eyepiece offers a certain magnification factor that has a determinative effect on the total magnification. Furthermore, all Leica eyepieces can be equipped with a direct image or a field of view measuring eyepiece (only on some models).

Diaptric connection

A diaphragm connection is also available for eyepieces when a more convenient viewing is required (see page 77).

If you do not wear glasses:

1. Hold the eyepiece firmly and rotate the eyepiece carrier towards the eye, as shown in the illustration.

2. If an eyepiece is equipped with the tilted diaphragm connection, turn the wheel in the direction shown.



If you wear glasses:

1. Hold the eyepiece firmly and rotate the eyepiece carrier towards the eye, as shown in the illustration. If the distance between the eyepiece and the eye is not sufficient, the eyepiece carrier should be moved further towards the eye.



2. If an eyepiece is equipped with the tilted diaphragm connection, turn the wheel in the direction shown.

By the way, one benefit of viewing with eyepieces is a distance of approximately 15 cm between the eye and the eyepiece (see page 76). The so-called eye relief distance ensures that your eyes will be relaxed even if they contact the eyepiece.



Focusing

Focusing is best achieved when the microscope is in coarse focus. The focusing drive mechanism should be engaged as the focus wheel is turned. The coarse focus wheel should be used for the initial focusing.

The focusing drive can be operated either manually or automatically.



Focusing

The fine focus wheel should be used for covering glass slides.

The outer, fine adjustment is used for fine focusing.



Coarse/fine adjustment

The sharpness is adjusted using the coarse/fine adjustment. The resolution of the coarse/fine adjustment is 1 μ m.



The coarse/fine adjustment carries a load of up to 15 kg.



Adjusting the Resistance of the Focus Drive

Adjusting the resistance

▶ The focus movement too loose or too tight?
Does the outfit tend to slide downwards? The resistance can be adjusted individually depending on the equipment weight and personal preferences to do so:

1. On the outer drive knobs with thumb wheels and turn their knobs clockwise until the desired resistance is reached during focusing.



Changing Magnification (Zoom)

4. After the camera level is leveled, zooming about the center indicates the zoom range covered:

- Zeiss M165 = 16.5:1
- Zeiss M205 = 20.5:1

The only camera parameter that needs to be set for zooming is:

- Zooming
1. Lock the keypad axes
 2. Focus on the subject first
 3. Rotate the magnification changer until the desired magnification is configured



Ratchet Steps and Magnification Levels

The zoom button can optionally be disabled. The user will not be able to adjust the ratchet steps or the zoom steps when the ratchet steps are disabled, which many users find convenient. On the other hand, when the ratchet steps are enabled, analogously, the zoom level buttons can be reproduced more accurately.

Enabling and disabling ratchet steps

1. Push the left button downwards to disable the ratchet steps.
2. Push the right button upwards to enable the ratchet steps.



Magnifications and fields of view
The user can [page 136](#) provides additional information about the magnifications and field of view details. It contains details on given to the position of the magnification changer and the eyes etc and object we can see when used.



Parfocality: More Comfort and Convenience for Your Work

- **Auto sleep** increases eye size by 10% when you close your eyes. When you open your eyes, the system automatically returns the eye size to the original magnification without having to refocus. The auto sleep needs only be re-enabled if you want to see a sheet that location that is located higher or lower.

Requirements for parfocal work

- If you're adjusting eye size and doing a correction, the procedure differs from this document for more information, refer to [page 73](#).
- For the procedure for adjusting the size of the object, the procedure refers to [page 29](#).
- For the procedure for adjusting the size of the object, the procedure refers to [page 33](#).

Parfocality

1. The eye never returns to the original eye size.
2. Focus on the sheet that

You're doing. Even if you see the sheet that you're adjusting the sheet that you're adjusting.

The size of the sheet that you're adjusting is the same as the size of the sheet that you're adjusting.



Iris Diaphragm

The iris diaphragm in the optics center of your M series microscope fulfills the same purpose as the iris diaphragm in cameras: It regulates the eye's depth of field, which changes the depth of field, the "depth of field" (or "focus depth") of the eye. It's a bit like the iris in your eye. It's OK.

Closing the iris diaphragm

- Close the iris diaphragm by turning the knob to the left. The subject appears darker and the depth of field increases.



Opening the iris diaphragm

- Open the iris diaphragm by turning the knob to the right. The subject now appears brighter, but the depth of field decreases.



Eyepieces



Magnification Factors of the Eyepieces

All eyepieces in the M series are designed to give the best view through the telescope, but also has a critical effect on the maximum magnification. The magnification factor is between 10x and 40x.

The following eyepieces are available for the M series:

Magnification	Optical correction	Order number
10x	1.5 d outer bellings	10 / 50 630
16x	1.5 d outer bellings	10 / 50 621
25x	1.5 d outer bellings	10 / 50 632
40x	1.5 d outer bellings	10 / 50 633



Health Notes

Potential sources of infection

⚠️ Direct contact with eyes, eyes, or
salivary secretions on the hand or
secretions in the mouth of the eye. The
secretions can be transmitted by using
eyes or direct contact with eyes. They can
be avoided by wearing eye protection.
Please contact your
doctor if they:



Separate eyecups are an effective
way of preventing infections.



Dioptric Correction

1. Set eyepieces and oculars to their full dioptric correction, ensuring the microscope to be used with your glasses even if those with vision problems. The correction corresponds to diopter settings.



Using the dioptric correction

1. Set the dioptric correction of both eyepieces to the red position (0 diopter) and the oculars to the blue position.
2. While wearing your glasses, look through the eyepieces and focus on the slide with the red arrow.
3. Rotate both eyepieces to the red position (0 diopter).
4. Hold one eye closed and rotate the other eyepiece in the red direction until the slide with the red arrow is visible.
5. Then, open the other eye and correct the diopter settings until the red arrow is clearly visible.



While the microscope is in diopter correction, the red arrow is visible through the eyepieces. To change the diopter correction, refer to the instructions on [page 78](#).



Dioptric Correction and Parfocality

Set the microscope eye piece to infinity. The microscope will be parfocal. The distance between the objective and the eyepiece will be the same. The distance between the objective and the eyepiece will be the same. The distance between the objective and the eyepiece will be the same.

Preparations

- Move the eye of the microscope to the distance between the objective and the eyepiece.
- You are adjusting the microscope to be parfocal.

Adjusting

1. Set the dioptric correction of both eyes to 0.
2. Select the lowest magnification and focus on a flat specimen.
3. Select the highest magnification and readjust the microscope.
4. Select the lowest magnification again, but do not adjust the eyes.
5. Take the eyepiece carrier out of the microscope and set it to infinity.

7. Look into the eyes.

8. Set the dioptric correction of each eye to the distance between the objective and the eyepiece.
9. Select the highest magnification and readjust the microscope.

Now if you adjust the magnification from the lowest to the highest, the subject will always stay in focus. This is the process.

Graticules

Use Leica graticules for length measurements and counting cells. Accuracy for measurements that are not equidistant: 1 in 10 µm.

The Leica graticules for length measurements and numbering are fitted in mounts and are inserted into the eyes.

1. Screw the insert onto the eyepiece.



2. Center the graticule on the insert, applying moderate pressure. Ensure that the graticule fits tightly.



3. Screw the insert and graticule firmly into the eyepiece.



4. You can move the graticule by rotating the eyepiece in the clockwise direction. Lightly rotating the eyepiece is recommended.




Use with the Leica AX carrier.



Use the Leica AX carrier with the Leica microscope. The Leica AX carrier is not suitable for other microscope models. The Leica AX carrier is not suitable for other microscope models.



Photography & Video



Photography & Video

For most measurement users, digital documentation is becoming the standard of the way. Test results can be presented in an effective manner; measurements on the digital image provide clarity and confidence. With the integrated software interface, even images of eye-bore holes can be captured and displayed on a computer for easier interpretation.

Adapter

Leica controls using the Leica AutoStation Suite is not required, conventional mirror reflex and rangefinder cameras from third-party manufacturers can be used. For this purpose, Leica Microsystems offers a variety of adapters that can be used together with the 30 A and 100 A Universal Units.

Leica DFC cameras

However, if you require accurate control over the camera and need the speed, flexibility and accuracy, the digital Leica DFC cameras are exactly what you need. Together with the Leica AutoStation Suite, they provide virtually all the features of a Leica camera, refer to the camera's documentation.



Leica Application Suite

The Leica AutoStation Suite is used to start a job. Here, the digital selection of the Leica AutoStation cameras, the addition of measuring targets, helps you control the measurement. In addition, images, camera and time for additional information, refer to the AutoStation Suite



Photo Tubes and C-mounts

Application

A C-mount camera lens equipped with a standard C-mount interface (1.27") and C-mount adapter for the respective photo tube is connected to this interface adapter. A cable and interface connection between the triaxial cable and camera ensures outstanding rendering of the triaxial image on the image sensor of the camera.

Due to the design of the digital camera to capture as much of the field of vision as possible, the video signal is truncated to the edge of the field of vision as possible. To do so, the magnification factor of the C-mount adapter must match the image angle of the sensor as close as possible.

The standardized video connection between the camera and C-mount adapter is made. It can be connected by the standard and/or of the camera software.

Alternatively, you can also use a C-mount adapter with higher magnification. This primarily avoids the critical border area of the field of vision and concentrates on the center of the field of vision.

Camera from third-party suppliers

In addition to the C-mount cameras with the standard C-mount interface, you can connect third-party cameras to the triaxial cable using a C-mount adapter.

To do so, instead of the C-mount adapter, you may use the corresponding S-1 adapter with a C-mount. However, these third-party cameras are not integrated into the camera software on Site and have to be operated using the corresponding software from the camera manufacturer.

The selected camera is deleted from the selected menu using the instructions for the camera's manual.

Time	Video sensor	Size (mm)	Field of view (mm)	Resolution (lines)	Gain	
9210040	5	29.76 x 19.23	2.37 mm	16	1633 40	0.5/0.5
9210040	10	29.76 x 19.23	1.07 mm	16	1633 40	0.5/0.5
9212001	3	29.76 x 19.23	3.2 mm	16	1633	0.5
9212001	3	29.76 x 19.23	3.2 mm	16	1633	0.5
9212007	2.8	29.76 x 19.23	3.5 mm	16	1633	0.6
9212001	5	29.76 x 19.23	3.5 mm	16	1633	0.6
9212001	5	29.76 x 19.23	3.5 mm	16	1633	0.5
9212005	8	29.76 x 19.23	3.7 mm	16	1633	0.6
9212001	12.5	29.76 x 19.23	4.25 mm	16	1633	0.6
92130015	1.3	29.76 x 19.23	3.75 mm	16	1633	0.6/0.6
92130015	1.3	29.76 x 19.23	4.25 mm	16	1633	0.6

© 2016, 2017, 2018, 2019, 2020, 2021, 2022, 2023, 2024

2025, 2026, 2027, 2028



Trinocular Video/Phototube 50 %

Use

With its third beam path, the trinocular video camera tube 50 % enables you to use the microscope video and to make use of the special third objective and eyepieces to view.

- 50 % beam path for the two eyepieces.
- 50 % of the light is diverted to the video camera beam path.

Assembly

Tighten the trinocular tube 50 % to the optics carrier block of the microscope observation tube. Refer to [page 21](#).



Trinocular Video/Phototube 100 %

Use

When using either the "trocular video" or "trocular 100 %" either you use the video or photograph the specimen. This means that 100% of the light goes either to the other side of the trinocular beam splitter or to the video deck.

Assembly

Fasten the "trocular tube 100 %" to the optics center block of the microscope observation tube. Refer to [page 221](#).



Switchover

- Turn the camera on the grills de of the tube into the horizontal position in order to guide the eye side grill into the eyepieces. You can now observe the specimen.
- Turn the camera on the grills de of the tube into the vertical position in order to guide the eye side grill into the camera. You can now photograph the specimen.



Trinocular Video/Phototube 100 %: ErgoTube 5° – 45°

Use

The Trinocular ErgoTube can provide the user a comfortable and ergonomic viewing position by allowing the user to fine-tune the viewing angle of the eyepieces to best fit their posture.


The swivel design allows the right user both 100% of light to the eyes as well as 100% of light to the camera port. The swivel allows the right user both to observe if they are not entering the tube.

Assembly

Before mounting the ErgoTube, please remove the Perspex block from under the tube. Otherwise the swivel design will be held.



Microscope Carrier



The Objective Nosepiece

Use

The objective noses extends as you rotate between two objectives just one hand movement, or else the objective is extended automatically.



Parfocal work

With the new Micro objectives, the distance between two objectives is automatically changed, meaning that the distance between the objective and the slide is always adjusting the height.



Order set objectives cannot be used, but you can take by using the objective height.

Technical constraints



The objective noses are cannot be used together in the 205000 Micro system. In addition, the objective noses are cannot be used.



Objectives and Optical Accessories



The Different Types of Objectives

To meet the various requirements regarding resolving power, there is a wide offering of fully achromatic, semiachromatic and semiachromatic objectives and also uncorrected, noncorrected objectives.

- Achromatic objectives are particularly suited for specimens with high-contrast structures.

- Hat-field (planachromatic) objectives are particularly well suited for studying flat objects such as cells and thin sections.

- With semiachromatic objectives, the finest structures are visible with high contrast. The so-called semiachromatic correction errors in these objectives result in the highest color brilliance and fidelity.

Achromatic objectives

The 0.32X, 0.5X, 0.63X, 0.8X, 1X, 1.5X, 2X achromatic objectives offer countless variants for selecting the object field diameter, magnification ranges and working distances (see page 156).


Planachromatic objective 1X

To the highest requirements for wide field of view, the semiachromatic objective with the 1X plan (flat-field) objective, which returns sharp, contrast-rich object fields,


Achromatic objectives with a long focal length

To observe easy colors, achromatic objectives with long working distances and low heights ($d = 100$ mm for 100 mm lenses etc.)

Bases




Leica TL ST Transmitted-light Base: Controls



1. *Adjusts plate to assembly/focusing dials*
2. *Removable glass plate*
3. *Control for light intensity*
4. *Adjustment for deflection mirror*



Extension plate for the transmitted-light Base TL ST



Back side of the transmitted-light Base TL ST

1. *Switch for changing the halogen lamp*
2. *Power connection socket*
3. *Power switch*



Leica TL ST Transmitted-light Base: Operation

Light intensity control

The LED control adjusts the intensity of the 12V/20W halogen light.

1. Select on the illumination of the base of the same section.
2. Turn on the base.
3. Set the illumination intensity using the LED control.



Transmitted-light control

The transmitted-light control automatically moves the deflection mirror of the base when moved. The transmitted light control also allows the deflection mirror to smoothly changeover between bright field and oblique transmitted light.



Bright field

Bright field is suitable for examining translucent objects featuring contrasting structures. The object is directly illuminated from above and is seen in its natural color against a black background.

- Move the deflection mirror. The desired effect is achieved.

Inclined transmitted light

The transmitted light illuminates the object obliquely. It is ideal for examining flat and thin material when observing very thin, translucent objects.

- Set the deflection mirror. The desired effect is achieved.



Leica TL ST Transmitted-light Base: Changing Bulbs

Changing the halogen lamp

Before you change the air, it is essentialy necessary to change the base for the base to prevent the base electric shock.

The halogen air becomes very hot during operation. Therefore, to avoid being burned, let the base cool off for approx 10 minutes after switching it off!

Do not touch the halogen air with your bare fingers – this drastically reduces the service life of the air.


Changing bulbs

1. Loosen the two screws on the feet sides and pull the feets out, along with the air.



Precautionary measures


When installing the air, ensure that the cables are inside the two cable covers. This prevents the cables from getting caught during operation.



2. Carefully pull out the air end through by pulling them upwards.
3. Disconnect the air from the mount.
4. Install the new air into the mount and re-install the air holder.




Leica TL BFDf Transmitted-light Base: Controls



1. Adjusts plate to easy-to-sensibility of focusing drives
2. Standard stage 10 (50 502)
3. Switch to toggle between bright field and dark field



Coarse adjustment of the transmitted-light base 71.9 mm



Adjusts of the focusing of fine




Connect to, for, cold light sources
 light control for, on line, 1 = 10 mm,
 and tube, 1 = 12 mm



Switch to toggle between bright field and dark field


Leica TL RC / Leica TL RCI: Controls



1. Heat sink of the integrated halogen illumination (only TL RCI)
2. Extension plate for focusing drives
3. Stator of stage 10 (50 502)
4. Filter holder
5. Control of top and bottom flaps of the filter mount controls
6. Button for fine and horizontal movement of the mirror
7. Transmitted-light base




Fig. 10a/cobase = Transmitted-light base
outlet for plate



1. Power switch
2. Power connection socket
3. IES socket, type 5
4. IES socket, type A
5. 2-CPM bus
6. Screws for changing the halogen lamp



Leica TL RCI: The Deflection Mirror



Tilted mirror

The built-in mirror features one flat and one concave side and can be tilted and moved. The concave side is intended to use at high magnifications to increase the NA of the base. The flat side is designed to deliver even an illumination at low magnifications. The black velocity knob on the left side of the rotary knob can be used to rotate the built-in mirror and move the field of view.

Inverted world?

Depending on the properties of the base, then the objective sides of the environment and the detection of the viewer, it is possible that the structures described here or will be and inverted real or virtual and to be detected in reverse. This means that the axes are not other than the upper section, camera inverted the real or virtual and is reverse.


Functions of the rotary knob

The rotary knob (1) fulfills the following tasks:

- Turning the deflection mirror from the flat to the concave side
- Setting it to guide the light beam through the base, then a circle sleeve or flatter angle
- Moving the deflection mirror (forwards or back)



Leica TL RCI: Color Intensity and Temperature



1. Potentiometer 1a controlling the **Color Intensity** Controls illumination intensity
2. Potentiometer 1b controlling the **Color Temperature**

The front bed-grip base **T** - **T** has two electronic potentiometers that control the color intensity, the color temperature (2).

The control of the color temperature is usually achieved by each electronic potentiometer:

- To alter just your work <<< the potentiometer (2).
- To change the position of the potentiometer, the electron returns the color temperature to the previously configured settings.

Using a USB mouse (only TL RCI)
The base USB mouse controls the **CC** end during the color **T** - **T** base. Connect the mouse to the corresponding USB port of the base.


- To alter the color of the images acquired by the TL RCI in the **CC** control system end is used to control the potentiometer by
- to switch the illumination on or off, briefly <<< the **CC** wheel.
- To enter or edit the **CC** end, click and hold down the **CC** wheel for longer than 2 seconds.



Leica TL RC / Leica TL RCI: Operation

Setting the relief contrast

The two switches on the left side of the built-in flap activate the relief contrast. The outer switch (1) controls the positive relief contrast, and the inner switch (2) activates the negative relief contrast.



1. Switch for adjusting the positive relief contrast
2. Switch for adjusting the negative relief contrast
3. Deflection mirror

Depending on the flap position, a part of the opening of the built-in lenses are covered, which results in the different contrast effects. Phase structures typically act as spatial relief contrast. The positive relief contrast sets the inverted relief contrast to the eye.



Increased contrast without relief is achieved by setting the deflection mirror to the illuminated area is created. By tilting the deflection mirror slightly, you can move the contrast over the entire field of view and quickly toggle between positive and negative relief contrast. The dynamic effect makes it easy to distinguish phase structures with the deflection mirror.



Leica TL RCI: Methods in Transmitted Light

Vertical bright-field illumination

Suitable for slide-mounted specimens with sufficient contrast

The light beam is deflected vertically through the specimen. This results in accurate bright field with maximum brightness

Inclined transmitted light

Suitable for semi-transparent, opaque specimens such as foraminifera and fish eggs. Move the deflection mirror until the desired data are visible

Singlesided dark field

Suitable for fixed specimens and fine structures

The flatter the angle at which the light beams are deflected into the specimen plane, the darker the substrate appears. A dark field-like transmitted light is created. Outlines, fine edges and structures are bright, in contrast with the dark background, through diffraction of the light beams on the dark background



Leica TL RCI: Relief Images

Starting position

1. Push the deflection mirror all the way back towards the camera.
2. Turn the deflection mirror into the notch position (see page 13).

Positive relief contrast

Slide to the left in the view and the view is dark. The structures are dark.

The effect can be strengthened or weakened by gently tilting the deflection mirror.



Negative relief contrast

Slide to the right in the view and the view is light. The structures are light.

The effect can be strengthened or weakened by gently tilting the deflection mirror.



Leica TL RCI: Relief Images (Continued)

Dynamic relief contrast

Subtle or very prominent and translucent structures

By tilting the deflection mirror slightly, you can move the gaps over the entire field of view and quickly toggle between positive and negative relief images. The dynamic effect makes it easy to distinguish several structures from each other at a glance.

Constraints

The relief methods provide good results for mid-zoom to high magnifications and with 1X, 1.6X and 2X objectives. The objective must be and with weaker specimens, the object field may not be uniformly illuminated.

We recommend using the Leica RCI-grip set with 1X or higher objectives and not objectives with a very low angle.




Using Filters

Filters for Leica TL RC and TL RC II

The Leica TL RC and TL RC II can be equipped with up to three filters – even at the same time. If a customer requests, the filters are also available as one-off items.

1. Switch off the light source or click the shutter button:



2. Take the empty filter from an available filter slot in the filter holder.

Daylight filter for Leica TL RC II
A daylight filter is also available for the Leica TL RC II
Leica Filter- grill case

3. Insert the desired filter.




4. Switch the light source back on.



Leica IsoPro (Non-motorized): Controls


Operating the Leica IsoPro mechanical stage

1. To move the stage in the X-Y direction, rotate the coarse knob.
2. To move the stage in the Z direction, rotate the fine control ring.




Leica IsoPro (Motorized): Controls

isoPro




1. Leica Soft or motorized mechanical stage
2. Housing with motorization

X,Y Stage DCI module



1. Sub-D interface to Leica Sma II/III
2. USB interface
3. USB Type B
4. Socket for grounded power cable


Leica PSC Controller



1. Quick control emergency function
2. Fine control in X direction
3. Fine control in Y direction



System Illumination



Leica LED5000 MCI

Use

The Leica LED5000 MCI for MJ11 Contact JTM theodolite is a universal right-angled JTM theodolite. Three groups of 3 LEDs each can be employed as desired around the base when the leveled control of the LEDs enters various JTM theodolite use is permitted.

Preparation

For the MCI with both the ends of 100mm rods, until the JTM theodolite is placed on the ball-bearing of the guide rods.

In this position, you always have the same contact with the contact JTM theodolite. This guarantees the vertical axis by a constant level.



Installed too high



Optimum height

Contact with the base



The optical centers are eccentrically oriented (see Fig. 2); contact between the base and MCI may result in safety-related or other irregularities. In this case, the MCI should be pushed upwards to prevent the base from being damaged.

- Be sure that the eccentric base is pushed fully beneath the MCI.
- After the optical centers' level, pull the MCI downwards near the base.



Leica LED5000 MCI (Continued)

Using the keyboard

- Use the **Q** key to lock the **ATT** function or off.
- Use the **+** and **-** keys to adjust the brightness in 10 increments. At the end of the list allows to adjust the intensity in five increments. Hold a key to change the intensity in a step.

The illumination scenarios


Use the **<** and **>** keys to toggle between the ATT illumination scenarios:

Maximum brightness

To increase the brightness, press the **+** key. The LED5000 MCI provides 10 increments of brightness. At the end of the list allows to adjust the intensity in five increments. Hold a key to change the intensity in a step.



The Leica LED5000 MCI can also be controlled by the LMS (Leica Application Suite) software.



Leica LED5000 RL

Use

The Leica LED5000 RL ring generates a very bright and homogeneous incident light. The diameter is 80 mm and the ring has a diameter of 20 mm. The Leica LED5000 RL can be used on and off completely or in various combinations. It can be used as a ring light or as a light source for the integrated LED.

LED enables you to create a very reproducible illumination scene. The Leica LED5000 RL can be used as a ring light or as a light source for the integrated LED.

Using the keyboard


- Use the **Q** key to switch the illumination on or off.
- Use the **+** and **-** keys to adjust the brightness in 10 increments. Touch either of the buttons to adjust the brightness in 10 increments. Hold the keys together to increase the brightness more quickly.

The illumination scenarios

Use the **←** and **→** keys to toggle between the illumination scenarios as shown below:


Constraint

The Leica LED5000 RL has been designed for use with the 1.0X and 0.63X objectives. For the 1.6X and 2.0X objectives, conventional ring holders from the Leica JTR holder product range are also used.



Leica EL6000 – About the Instrument

The Leica EL6000 compact gill source is intended for surgery and the end section works without the use of a sectioning table. The turner body of the source is getting smaller.

 The Leica EL6000 compact gill source is provided exclusively for operating mercury switches of type 402-1120W/25 VS from OS74M in conjunction with the source gill guide for the instrument. Any other use of the instrument is to be considered non-conformant use.

Special features

The Leica EL6000 compact gill source has an automatic beam cover. This ensures that the gill beam outside of the instrument when the gill guide is detached. This prevents the beam from being harmed by the intense gill beam generated by the instrument.

The Leica EL6000 compact gill source has a built-in shutter. It can be operated at the instrument by means of a push button or via remote control. The shutter is automatically closed when the user is not needed.

The gillness can be varied in 5 increments (0%–100%) for this source. The reference driving unit can be operated using a rotary switch on the front of the instrument in the device.

Heat-absorbing filter

When the cover is open, the window of the heat filter can be accessed. Here, (hardened) heat-absorbing filters for protecting the light guide and various conversion filters with a diameter of 32 mm for adapting the beam from the excitation source are inserted.

The Leica EL6000 comes factory-equipped with a heat-absorbing filter.

Leica EL6000 – About the Instrument (Continued)

Compatible light guides

The Leica EL6000 camera grid source can be connected with grid light guides that are compatible with the Star and grid light end whose transmission curve fits the lamp used. Whenever possible refer to the manufacturer's or baseplate grid data sheet.

When connecting the grid guide to the grid source of the microscope edition, make sure that it is correctly fitted to the grid guide.

Your Leica camera lens is designed to provide additional information about baseplate grid guides.

Identification of the instrument

The reference of the instrument is on the bottom. It includes the necessary information for identifying the instrument.




The Leica EL6000 camera grid source can be connected with grid light guides that are compatible with the Star and grid light end whose transmission curve fits the lamp used. Whenever possible refer to the manufacturer's or baseplate grid data sheet.



Leica EL6000 – Operation

Controls



1. Power switch/light on indicator
2. Reset button for running-time meter
3. Running-time meter
4. Stop the control
5. Intensity switch

Running-time meter

The running-time meter measures the running time of the circuit that the system from the Leica EL6000 control unit source. The running-time meter displays a current value in this circuit hours - minutes.

1. Press the reset button to set the running-time meter of the circuit to '0'.

The running-time meter should be reset after every circuit measurement so that it always displays the current, accurate circuit value.

3 scroll the running-time meter when it reaches the end of the life specified by the manufacturer.

Shutter control

The 'user' position the shutter of the Leica EL6000 is open – regardless of the signal from the remote control.

The shutter control is the 'remote' position, the shutter is controlled by the remote control.

Intensity switch

The intensity switch provides a fixed amount of light in five fixed increments.



Leica EL6000 – Troubleshooting

Fault description	Possible faults	Remedy
Power indicator signal not indicated after being switched on, instrument shows no function	No power supply voltage present; connecting cable not plugged in correctly/defective - defective instrument defective	Check power supply voltage and connecting cable Test cable, use certified replacement If necessary, see page 114
Power switch is not indicated after being switched on, instrument shows no function	- faulty power cable - power switch - Stuffer closed or not fully inserted in slot	- Check for correct connection - Test cable - Open stuffer or rotate into position - Insert into slot
Target system receives too little light or no light	- grid guide not inserted into grid guide correctly - grid guide damaged/defective	- Check grid guide - Test or replace

If these actions do not fix the problem, send the Leica EL6000 compact grid guide to Leica Service for the exact description of the malfunction.



Cleaning and Maintenance



Be sure to carry out cleaning tests, switch off the Leica E6000 compact grid source.

Be very edify. It takes a long time when cleaning. Ensure that you do not alter the grid and enter the device through the ventilation slots.

In case of damage, even extreme damage, immediately stop using the Leica E6000 compact grid source and send it to service.

The manual should be consulted for the Leica E6000 compact grid source for details about the safety cover, such as the safety bag for the grid source setting.

Spare parts and accessories

The only parts of the Leica E6000 compact grid source that may be replaced by the user are the instrument cables and the eras.

- The replacement of the eras is described on [page 62](#).

Replacing the instrument fuse

1. Disconnect the Leica E6000 from the power supply.
2. Open the device in the instrument case described in the user's safety manual.
3. Test the instrument case.

Permitted replacement fuses

-K100/20, 2.0 A, 5.0 V, 5.0 W, 10 A, 250 V, 1000 mA

Manufacturer's: www.vishay.com

Supplier's: www.schurter.com

Repairs

Opening or working on the Leica E6000 compact grid source may render other than described above (changing eras, changing fuses) that may be carried out by the user or the adjustment of the source or other adjustments themselves.

Fluorescence





Filter Changer

The patented FLUOVIEW filter system consists of a filter changer for barrier and excitation filters, a 4-position filter carousel and a filter holder with a choice of filter.

The filter changer can hold a total of 4 filter sets. The filter sets for fluorescence are labeled (see page 119). The fluorescence filters have a transponder with the specific filter data for the Leica M165 TC, M205 TC and M205 TC.

The size of the Leica M205 TC, M205 TC and Leica M165 TC includes three simple filter holders without fluorescence for empty filter positions and transmitted or incident light observation without fluorescence. These filter holders have two empty openings for the barrier filter sets. The filter sets are labeled (see page 119).

 If you are using fewer than four filter holders with fluorescence filters, always insert the included simple filter holders with closed barrier filter sets into the carousel. Otherwise, there is a risk of damage to the eyes from direct illumination from the barrier filter sets.

 The simple filter holders also protect the system from dust.

Filter sets for fluorescence
The filter sets for fluorescence contain, on one filter holder, two barrier filters to the visual beam paths and one excitation filter to the illumination beam path. The fluorescence filters have a transponder with the specific filter data for the Leica M205 TC and M205 TC. The filter sets are labeled (see page 119).



A selection of the available filters including simple filter holder (bottom right)

About Fluorescence Microscopy (Only Applicable for Leica M205 FA)

Functional principle

Fluorescent substances emit light when excited with short-wave excitation light, and fluorescence microscopy takes advantage of this property. Specimens without self-fluorescence are stained with a fluorescent substance. One such dye is green fluorescent protein (GFP), which is used in molecular biology.

The Leica M Series

The Leica M series fluorescence stereo microscopes allow entire fluorescent specimens to be studied, unexcised and in three dimensions. The fully automatic auto system and the ergonomic work range make the Leica M series stereo microscopes your first choice. The patented FLUODIP filter system offers fluorescence illumination and the patented FLUODIP filter system generate highest-quality fluorescence images. The very bright high-pressure mercury burner and specially selected filter sets enable you to distinguish the finest structures and enhance the quality of illumination for incident light fluorescence.

The filter system

The FLUODIP filter system consists of a rapid filter changer for barrier and excitation filters and a filter insert for a filter holder with your choice of filter. A total of four filter sets (which are labeled) can be used in the rapid filter changer. An empty filter cartridge is also available for individual filter combinations. Each filter set includes two barrier filters to the visual beam paths and an excitation filter to the illumination beam path.



Rapid Filter Changers and Filter Types

The rapid filter changer must be equipped with our barrier filter models etc. Please refer to the standard filter sets table for details.

Filter sets	Excitation filter	Barrier filter	Designation
G ⁰¹	230/0 1R 505-550 1R1	430_2	G ⁰¹
G ⁰² JS	430/0 1R 760-800 1R1	510_2	G ⁰²
G ⁰³ JETS	470/0 1R 750-800 1R1	525/50 1R 500-550 1R1	G ⁰³
JV	560/0 1R 570-580 1R1	720_2	JV
Voet	230/0 1R 705-775 1R1	660_2	V
3 JE	470/0 1R 750-800 1R1	515_2	3
Gree 1	570/0 1R 571-581 1R1	590_2	G
C ⁰¹	560/20 1R 720-770 1R1	430/0 1R 760-800 1R1	C ⁰¹
Y ⁰¹	510/20 1R 500-520 1R1	560/0 1R 570-580 1R1	Y ⁰¹
Texas Red	560/0 1R 570-580 1R1	610_2	TR
DsRED	575/50 1R 510-560 1R1	620/60 1R 590-650 1R1	DS
C15	620/60 1R 590-650 1R1	400/5 1R 360-450 1R1	C15
C15	555/60 1R 550-580 1R1	610/5 1R 575-613 1R1	C15



Simple Filter Holders

Simple filter holders

Three simple filter holders are included with the Jetco M165 T and the Jetco M205 T. The M205 K.A. these filter holders have two empty openings for the observation beam light and opening to the ultraviolet beam light as covered.

- Use this filter holder if you want to work briefly (no more than 15 seconds) without fluorescence illumination.



So that the filter holder is not heated, ultraviolet does the mercury grid for longer than 15 seconds.

Safety Notes

- If you want to look the mercury grid or angle the filter holder the Jetco filter holder.
- Use the filter if you want to change to transmitted grid or to direct grid observation without fluorescence.
- Always place a simple filter holder into the empty filter changer positions in order to protect yourself from ultraviolet and scattered light.



If you are using fewer than four filter holders with fluorescence filters, always insert the included simple filter holders with covered ultraviolet beam light into the eye protection. Otherwise, there is a risk of danger to the eyes from direct ultraviolet from the UV beam light.



Equipping the Filter Changer



When inserting the filter sets, avoid touching the filters in order to avoid fingerprints. Clean dirty filters immediately by soap, water and alcohol.

Removing the simple filter holder

The Jetec M165 and the Jetec M205 and M205 T4 are equipped at the factory with simple filter holders as "dummies". These must be removed before inserting a filter.

1. Grasp the recess of the dummy in the filter changer and pull out.



Inserting the filter

1. Place the filter into the filter changer.

Ensure that the contour of the filter holder matches the contour of the recess.

2. Rotate the filter changer by hand until it clicks into place.

If the filter holder has not been inserted correctly, the filter changer does not rotate.


3. Insert a total of four filters.


Slot for filter slide


There is a slot in the filter changer system for an individually selectable filter, for example a neutral density filter.



Observation Without Fluorescence

 We recommend that you first familiarize yourself with the controls of your wide-field fluorescence microscope in transmitted light without the presence of fluorescence illumination.


 After switching on the high-pressure mercury burner, you must adjust the excitation filter wheel (see Section 12-216-11).

 We recommend that you first familiarize yourself with the basic functions of your wide-field fluorescence microscope before the observation of the specimen. These include the adjustment of the eyepieces, the zoom magnification changer, the focusing drive and others.

The following describes the specific functions of the fluorescence system only.

Observation without Fluorescence

The FLUOfilter system includes a UV shutter that can be used to observe the specimen without the UV light. This feature protects the mercury emission, which should not be switched on and off too frequently. In addition, you need to interrupt your work; you do not have to switch off the high-pressure mercury burner to observe the specimen. The UV light is being turned off automatically.

 When working with the 26000, it is recommended to use the shutter on the emission side of the grid guide. Switching the high-pressure mercury burner on and off frequently decreases its life and causes deposits on the emission side. Switched back one later than usual. Instead, you can use the UV shutter and reason 1 when you return to your work etc.

Close the UV shutter when you do not want to view the specimen with fluorescence illumination. Autogen: H11151711ed - gnt.



Double Iris Aperture

The built-in double iris diaphragm regulates the depth of field. Note that

- the greater the depth of field, the darker the image becomes.
- if the same size, the exposure time becomes longer or analogically.
- The resolution decreases.





Adjusting the double iris diaphragm using the knob below.



Commissioning the Fluorescence System

After switching on the high-pressure mercury source, you must adjust the mercury intensity. See [Mercury Menu \(M2-216-1\)](#). Do not switch on the supply unit until you feel confident using the adjustable mercury source.

 To avoid safety hazards when handling mercury cells and the cells housing

 For observation with fluorescence light, always position the probe on a screen correctly to protect the observer. See [page 271](#).

1. Set the mercury intensity. See [Mercury Menu \(M2-216-1\)](#).
2. Wait 2-3 minutes.
3. Adjust the mercury recorder with [Mercury Menu \(M2-216-1\)](#).



Accessories



Leica Hand Wheel and Foot Switch



Up to five foot switches can be configured with the microscope system. They can be programmed in LAS or LAS X to control focus, filter wheel, Z-position, and individual Z-steps.



The hand wheel is designed to control the motor.



Leica SmartMove



With its smart controls and smart, the Leica SmartMove control unit offers the right tool for handling motorized adjustments.

Controls

The functions of the individual controls can be configured both directly on the Leica SmartMove and also by the Leica Application Suite.

To see detailed description of the functions and settings, please refer to the Leica SmartMove manual or the Leica Application Suite help file.



Leica SmartTouch



With the touch screen of the Leica SmartTouch, you can control your edge point traverses with the touch movements and never use a grid or your mobile distance meter.

The most important control functions on the control unit can be adjusted to your specific needs with the key programmed safety status and junction buttons.

Controls

The functions of the individual status and screen elements can be configured both directly on the Leica SmartTouch and using the Leica Assistant Suite.

The standard configuration when shipped from the factory is shown:

- Top left of the status = count
- Top center of the status = focus
- Buttons = left filter, right filter, shutter

To delete description of the junction and settings, please refer to the Leica SmartTouch manual or the Leica Application Suite help file.



Leica PSC Controller

Use

The Leica PSC controller gives you precision control of the Leica SP16 motor and rectilinear stage.

Connection

Connect the Leica PSC controller and the Leica SP16 motor and rectilinear stage to the DX1 box.




Passing over the specimen quickly

- Move the optics in any desired direction, adjusting diagonally to control the rectilinear stage quickly.



Control by LAS

The instructions for controlling the Leica SP16 by the Leica Application Suite are provided in the software's help file.



Fine control of the mechanical stage

The motorized Leica SP16 rectilinear stage offers an accuracy of up to 0.25 µm. To move to a sub-micron level accuracy, the fine control stage of the optics




Dimensional Drawings




Leica M165 FC

Leica M165 FC with transmitted-light stand and focusing column
(dimensions in mm)






Leica M165 FC (Continued)

Leica M165 FC with transmitted-light stand TL RC, manual cross-stage, Leica bu-Pro and manual focus




Leica M165 FC (Continued)

Leica M165 FC with transmitted-light stand TL RC and manual focus




Leica M205 FA

Leica M205 FA with transmitted-light stand TL BFD F, trinocular EyepTube, DFC camera and motorized focus
(dimensions in mm)




Leica M205 FA (Continued)

Leica M205 FA with transmitted-light stand TL RD, binocular ErgoTube, DFC camera and motorized focus






Leica M205 FA (Continued)

Leica M205 FA with transmitted-light stand TL RD, IsoPro mechanical stage, trinocular ErgoTube, DFC camera and motorized focus




Leica M205 FCA

Leica M205 FCA with transmitted-light base and focusing column
(dimensions in mm)






Leica M205 FCA (Continued)

Leica M205 FCA with transmitted-light stand TL RC1, Leica body manual mechanical stage and manual focus






Leica M205 FCA (Continued)

Leica M205 FCA with transmitted-light stand TL RCI and manual focus





Leica LED5000 MCI

Leica LED5000 MCI
(dimensioni in mm)





Leica LED5000 RL

Leica LED5000 RL
(dimensionen in mm)



Leica MATSTL



Specifications



Leica M165 FC / Leica M205 FA / Leica M205 FCA

	Leica M165 FC	Leica M205 FA	Leica M205 FCA
Optical data			
Zoom	16.5:1 variable, coded	20.5:1 variable with zoom dial	20.5:1 variable, coded with zoom dial
Delivery intended optics (1X object view/10X eyes view) - Zoom range - Resolution - Working distance - Object field	7.2X-120X Tex./52 µm 61.5 mm, 3 telephoto lens Ø 31.5 mm - 192 mm	7.8X-160X Tex./525 µm 61.5 mm, 3 telephoto lens Ø 29.5 mm - 177 mm	7.8X-160X Tex./525 µm 61.5 mm, 3 telephoto lens Ø 29.5 mm - 177 mm
Viewing device, based on optics construction - Magnification - Resolution - Visual acuity width - Entrance aperture - Object field	960X 906 µm 551 mm 0.302 Ø 63 mm	1,280X 1,050 µm 76 mm 0.35 Ø 59 mm	1,280X 1,050 µm 76 mm 0.35 Ø 59 mm
Working distances		135 mm, 0.5X telephoto lens 112 mm, 0.8X telephoto lens 67 mm, 0.63X telephoto lens 61.5 mm, 1X telephoto lens 30.5 mm, 1.6X telephoto lens 20.1 mm, 2X telephoto lens	

Leica M165 FC / Leica M205 FA / Leica M205 FCA (Continued)

	Leica M165 FC	Leica M205 FA	Leica M205 FCA
Optics carrier			
100% electronic distance system	CMO Common Metric Objectives only		
Specific surface resistivity (microhm)	2x1,011 Ω/mm^2		
Coded/decoded; M205 FA	distance line (2 seconds) from 1,000 V to 100 V		
Coded/decoded; M205 FA	Zoom, 1/3rd integr., objective masks etc.	Zoom, 1/3rd integr., objective nosepiece, FM filter	Zoom, 1/3rd integr., objective nosepiece, filter
Cybercode zoom markings	1/3rd reset live lock	1/3rd reset live lock	1/3rd reset live lock
Distance 1/3rd integr. for detail of field control	3rd 1/3rd coded	3rd 1/3rd manual	3rd 1/3rd coded
1/3rd 1/3rd reference intensity control	-	3rd 1/3rd manual	-
Accessories			
Single objective take off the objective masks etc.	1x 3rd reset manual		
Full time objective take off the objective masks etc.	5x 3rd reset manual, 2x 3rd reset manual, 1.6x 3rd reset manual, 0.62x 3rd reset manual		



Leica M165 FC / Leica M205 FA / Leica M205 FCA (Continued)

	Leica M165 FC	Leica M205 FA	Leica M205 FCA
Additional objectives, wide-field eye	1X objective, 0.8X objective, 0.63X objective, 0.5X objective, parallel		
Interchangeable eyepieces	50 mm – 105 mm		
Wide-field eyepieces for persons wearing glasses	10X, 16X, 25X, 40X, 71 mm eye relief, 3.5 & 4.5 diopters		
Manual coarse/fine focus	Focus range 130 mm or 330 mm, adjustable movement	Focus range 130 mm, adjustable movement	Focus range 130 mm, adjustable movement
Motor and focus drive	Motor and focus drive 11/21 mm profile column, with power supply		
Computer interface	J63		



Leica M165 FC / Leica M205 FA / Leica M205 FCA (Continued)

	Leica M165 FC	Leica M205 FA / Leica M205 FCA
Module System	Stands, illuminators	
Transmitted light stands	For bright field and bright/dark field, high-performance base: Bright field, single-sided dark field, with Turret in Control and CCD	
Reflected light stand	Large x-height grid stand with accessories like stage elevators and shutters	
Stages	Grid stage, MTS feeding stage, carriage, table elevator stage, stage support, microscope stage, turret and motor feed	
Light sources	LED 5000 MCI and LED 5000 RL, fluorescence, oblique, coaxial, cold light sources	
	Misc. accessories	
Photography video	<ul style="list-style-type: none"> Leica Digidig microscope system, microscope video Video adapters for camera, video and CCD camera 	
Image capturing, image processing	Leica AFluor: Special software for all fluorescence tasks Leica AFluor, Leica AFluor, Leica AFluor, Leica AFluor, Leica AFluor, Leica AFluor	
Measurement grid cues	For high-resolution counting	
Vertical and/or horizontal	Z-axis detector and microscope elevators	
Directing laser	For both reflected and transmitted	

To delete information about the accessories, please refer to the Leica M205 FA or M205 FCA manual. For more information, see the Leica M205 FA or M205 FCA manual. For more information, see the Leica M205 FA or M205 FCA manual.



Leica EL6000 – Specifications

Dimensions (WxHxD)	
Width x height x D	120 mm x 210 mm x 290 mm
Environmental conditions	
Permitted ambient temperature	-20 to 85 °C
Permitted relative humidity	10% to 90%, non-condensing
Operation	
Permitted ambient temperature	0 to 40 °C
Permitted relative humidity	10% to 90%, non-condensing
Operational data	
Accessories	Discussed separately
Warranty for use	Max. 2,000 measurement cycles
Protective class	
Type of protection	IP20 (EN 60529)
EMC	
Electromagnetic interference immunity to the device	EN 61326 Class 31
Voltage fluctuations, flicker	EN 61326, Annex A4
Electromagnetic compatibility	EN 61000-3-3
	EN 61000-3-2
The instrument fulfills the requirements of EC Directive 89/336/EEC and the EMC standard September 18, 1998	

Electrical safety	
Standard reference: IEC 61010-1:2001	
Protection degree	2
Overvoltage category	
Power input	Selected socket: EN 620 C13 EN 60320 C13
Design	
Rated supply voltage	100 to 240 VAC, 50/60 Hz
Rated frequency	50 to 60 Hz
Tube protection	
Tube protection	2x 2.5 A, 3 times triggered by 4. The manufacturer: Wüchtrich Series 1811
Power consumption	Max. 210 VA
Lamp	
Type	Mercury short arc reflector lamp
Designation	OSRAM-HRS-T120W/2EVS
Power consumption	120 W
Service life (manufacturer's information)	Average life 2,000 h



Leica EL6000 – Specifications (Continued)

Light output	
Design	Type Star, ring
Socket feature	Set locking into gridguide is jogged in
Fan	
Switch-on delay	Absmax. 60s after plugging the lamp
Interface	
Type	Terrate control terminal
Design	9-pin D-sub socket
Terminal type	
- Stutter closed	-8 to +0.8 V
- Stutter open	+2.7 to +8 V
- Fault current	Absmax. 1 mA
Footswitch	
- Stutter closed	Continuation between 5 and 7
- Current	Absmax. 1.2 mA
Stutter frequency	Max. 10 Hz

Protective functions	
Safety cut	In the case of a short circuit, the lamp current is cut off. In the case of a fault and power is again connected, the lamp supply voltage.
Overheating protection	In the event of overheating, the air supply is shut off. After cooling off, the lamp power supply in the case EL6000 cannot be restarted automatically.
Protection screen	In the gridguide is jogged into the gridguide terminal, the gridguide seal can be closed off so that no light can escape.



Leica TL ST Transmitted-light Base

Light source	Halogen lamps 12V/20W
Quick illumination change	Yes
Unfiltered light	50 nm
Power supply unit	Autova Light 100–230 V, frequency 50–60 Hz, energy consumption 30 W max, ambient temperature 10–30 °C
Connectors	Power supply, power switch
Weight	7.7 kg
Illumination types	
Bright field	Yes
Dark field	Yes (integrated)
Oblique light	No
Leica Contrast System (LCS)	No
LCS Contrast System (LCS II)	No
Mercury filter/attachment	No
Integrated filter holder	Yes
Coated objective increasing the contrast	Yes
Temperature	No
Making a ring illumination	No
Temperature variations	No
AntiShocks Pads	Yes
Dimensions (WxHxD)	270x/30x85 mm



Leica TL BFDf Transmitted-light Base

Light source	Mercury cadmium
Lighted area	20 mm ²
Connectors	Connection card guide, cable = 10 mm, end tube = 13 mm
Weight	5.8 kg
Illumination types	
Bright field	Yes
Dark field	Yes
Oblique light	No
Leica Contrast System (LCS)	No
LCS Contrast Control Intensity Control	No
Mercury anti-reflecting coating	Yes**
Integrated filter holder	No
Coated objective lenses by the coating temperature	No
Making of microscope	Yes**
Temperature variations	Yes***
Anti-static pads	Yes
Dimensions (WxHxD)	370x290x90 mm
*With card guide source Leica CS150 S	**Contrast control
	***With mercury guide source



Leica TL RC / Leica TL RCI

	Leica TL RC	Leica TL RCI
Light source	Helium neon grid source	Helogen lamp 12V/20 W
Quick attachment	-	Yes
Interleaved lines	35 mm	35 mm
Power supply unit	-	100-270 V, frequency 50-60 Hz energy consumption 30 W max, ambient temperature 10-30°C
Connectors	Connection board, grid guide, active = 10 mm, end tube = 13 mm	1 x J53 type A, 1 x J53 type B, 2 x C.A. 4-3 J5
Weight	6.0 kg	7.2 kg

Illumination types

Bright field/dark field	Yes/yes	Yes/yes
Optical grid/Telet Control System (TCS)	Yes/yes	Yes/yes
TCS Constant Contrast Control	No	Yes
Telet Control/Contrast Control	Yes**	Yes
Integrated filter holder	Yes	Yes
Control of increasing the contrast	Yes	Yes
Masking of illumination	Yes***	Yes***
Telet Control actions	Yes****	Yes
Anti-Streak Beds	Yes	Yes
Dimensions of base, WX-HxD	370x390x95 mm	370x770x95 mm
*Singlesided **W in card, grid source, base C.S130 S ***Contrast ratio ****W in electric grid source		



Leica TL5000 Ergo

Light source	
Light source	LED
Included accessories: 3-year warranty	65 mm
Included accessories: Desk stand	40 mm
Teletest Contrast System (TCS)	Yes
Merle Staller's microscope	Yes
Teletest contrast system	Yes
Included letter reader	Yes
Making digital microscope	Yes
Print-Snap & Peds	Yes
Size (WxD) (mm)	412 x 341 x 46
Accessories	
Software: ImageSight Support	-
Included test rig: Single / dual / triple	Yes
Included test rig: double / triple / single	Yes
Ergo Test / dual / triple	-

Power supply	
Input	120-240 VAC 47-63 Hz 3.2 A
Output	33 VDC 3.95 A 100 W max
Connections	
COM	3
USB	1
Power Supply	1



Leica IsoPro Motorized XY Stage

Control type	Leica intelligent base and base
Stages size (X×Y×H)	335.5 mm×270 mm×15 mm
Travel	152 mm×102 mm×6 mm
Speed XY	30 mm/s
Motor resolution (X,Y)	0.25 μm
Resolution (X,Y)	1.2 μm
Positioning accuracy (X,Y)	±20 μm over the entire positioning range
Types (X,Y)	30 μm
Motor	Direct current (DC) motor with encoder
Accessories	See text
Indicators	Indicators for base
Acceleration	Systemed acceleration
Controller connection	USB
Control options	Leica PSC, Leica Steel Move, Leica Steel Touch, Leica iSive PC
Power supply	Input: 100–240 VAC 50/60 Hz 1.0 A Output: 15 VDC 2.66 A/10 W max.



Appendix



Calculating the Total Magnification and Field of View Diameter

Parameter	
MO	Magnification of the objective
ME	Magnification of the eyepiece
Z	Magnification of the changer position
q	Tube factor; eg 1.5X or case incident light 1.6X or 1.5X eye-tube
.	Tube factor 1.25X of the binoculars and the normal objectives of the M125/M116 are used on the M55, M26, M275 or M295
FTOV	Field number of the eyepiece. Field numbers are field on the eyes. eg 10X = 21, 16X = 17, 25X = 9.5, 40X = 6.

Parameter	
MO	1X objective
ME	25X/9.5-eyepiece
Z	Zoom position 1
q	Coaxial reflected light 1.5X tube factor
.	Factor 1.25X

Calculation example: Magnification in the binocular tube:

$$\text{MFTVS} = \text{MO} \times \text{ME} \times \text{Z} \times \text{q} \times .$$

$$1 \times 25 \times 1 \times 1.5 \times 1.25 = 187.5 \times$$

Calculation example: Field of view diameter in the specimen

$$\text{FOV} = \frac{\text{FTOV}}{\text{MO} \times \text{Z} \times \text{q} \times .}$$

Care, Maintenance, Contact Persons

We hope you enjoy using your grease gun. Grease guns, grease guns, grease guns are determined by the manufacturer and are not covered. Observing the operating and maintenance instructions will ensure that even after years and decades, your grease gun will continue to work as well as it did on the very first day.

Warranty benefits

The warranty covers electrical failures and malfunctions. It does not, however, cover damage resulting from misuse or improper handling.

Contact address

For more information, contact your distributor or, if you are unable to reach your distributor, contact your local office. Write to: info@greasegun.com, or call 1-800-451-1111.

E-mail contact

info@greasegun.com

Care

- Protect your grease gun from moisture, dirt, and debris. Store in a clean, dry place. Do not use the grease gun for anything other than its intended purpose.
- Do not use the grease gun for anything other than its intended purpose. Do not use the grease gun for anything other than its intended purpose.
- Do not use the grease gun for anything other than its intended purpose. Do not use the grease gun for anything other than its intended purpose.
- Do not use the grease gun for anything other than its intended purpose. Do not use the grease gun for anything other than its intended purpose.

Care, Maintenance, Contact Persons (Continued)

Protection from dirt	Cleaning polymer components	Permitted measures
<p>Dust and dirt will affect the quality of your results</p> <ul style="list-style-type: none"> • The dust cover over the microscope when it will be used during the • Use dust caps to protect tube openings. Use only dust caps, and eyes etc. • Feet accessories. The dust-free is also when not in use 	<p>Some components are made of polymer or are polycarbonate. They are therefore essential and convenient to avoid the use of abrasive cleaning agents and techniques that damage polymers.</p>	<ul style="list-style-type: none"> • Clean the microscope parts of all dry matter by water. The microscope dust adheres: • For adjustment, you can use alcohol and the ethanol for adjustment. When doing so, use the corresponding safety instructions. • The removal of the microscope is possible only by the manufacturer. • Clean the objectives and eyepieces with special dust cleaning cloth and in the camera.





CONNECT
WITH US!



Leica Microsystems Schweiz AG · Max-Born-Strasse 201 · 9135 Heerbrugg, Switzerland
+41 71 726 3131 · +41 71 726 3141

www.leica-microsystems.com