Living up to Life



Leica M720 OH5

A Paradigm Shift in Vision, Comfort, and Flexibility Premium Surgical Microscope for Microsurgery

A Paradigm Shift

- Comfort through ergonomic design
- Brilliant images
- Unsurpassed patient safety
- Intraoperative fluorescence
- Viewing for the entire OR team
- Positioning flexibility
- Superior maneuverability





Leica M720 OH5 A Paradigm Shift in Vision, Comfort, and Flexibility

For years, surgeons have needed a surgical microscope with smaller, more compact optics. Traditional microscope design has evolved over the years using large, vertical optical zoom lens systems, which have inherently limited the surgeon's amount of working room, and the ability to work in the right ergonomic position. With the Leica M720 0H5, Leica Microsystems writes a revolutionary new chapter in microscope innovation. At the heart of the innovation: Horizontal Optics Technology.

> The heart of the innovation: Horizontal Optics Technology reduces the size of the optical head and gives the user more space to work. At the same time it dramatically increases comfort.



More Space to Work

Leica M720 OH5 Comfort Through Ergonomic Design

The Leica M720 optical head is the most compact of all neurosurgical microscopes. Designed along a horizontal plane, the compact optics carrier helps the surgeon naturally align for a healthier working posture. Whatever the position of the patient, even sitting upright, the surgeon can see more, can work more precisely, and benefits from superior ergonomics.

Compact Horizontal Optics

The substantial gain in free working distance gives the surgeon unobstructed access to the surgical area, greater instrument maneuverability, and an optimal view.

Butterfly Binoculars


Leica's butterfly binoculars accommodate all body heights, for both the surgeon and the assistant, as well as the most challenging surgical positions. The tubes have an inclination range of 115°, and the eyepieces swing to a second viewing plane quickly and easily.

SpeedS po t[™]

Two laser beams act as a focusing reference to quickly provide a defined focus point for all three viewing ports (surgeon, assistant, and camera).










Light Where You Need It

Leica M720 OH5 Brilliant Images

The Leica M720 OH5 is equipped with Small Angle Illumination (SAI) to distribute more light to the bottom of deep cavities. SAI provides a concentrated light beam, closely aligned to the optical axis.

Combined with outstanding Leica APO OptiChrome[™] optics, the result is significantly improved depth perception and better light penetration, specifically for new surgical access techniques such as intra-tracheal, transsphenoidal or METRx[™]. Images have outstanding contrast, brilliance, sharpness, resolution, and color fidelity.




e.g. transsphenoidal surgery, general illustration

Enhanced 3D Images: Depth perception is improved thanks to Lore: elargo etoro ob aco of D4 mm creating a more true-to-life 3Deffect compared to other microscopes.









Leica M720 OH5 Unsurpassed Patient Safety

Dual independent light sources [1] : The Leica M7200H5 features two completely

inde pendent 400 W xe non a rc-l amp illumination systems, givin g the surge on confidence to know that surge ry will not be jeo pardized due to lamp or bo ard fai lure.


Fast system reboot: If the power cable be comes disconnected for any reason, the system reboots in the fastest reset time available to day.

In tuitive user controls [2]: The graphical user interface and hard keys allow the user to conveniently and intuitively control all microscope functions during surgery.

In dependent microscope controls: Stand, video, light, and microscope controls work independently. For example, should the video system fail, surgery can continue because the light and microscope are un affected.


Antimic robial surface coating [3]: Leica's Ag Protect™ limits user exposure to surface pathogens. This nano silver coating covers the miαroscope's outside surfaces and pen etrates the membranes of microbes to prevent replication.

The Leica M720 OH5 offers in novative illumination solutions to improve outcomes for both the surgeon and the patient.





Conven tional



BrightCare Plus – Light Intensity

BrightCare Plus optimizes the light intensity relative to the working distance. As working distance decreases, the light intensity is reduced automatically, minimizing incidents of patient burns. As working distance increases, the light intensity rises again accordingly.

Autolr is[™] – Light Diameter

Auto Iris[™] automatically adjusts the diaphragm so that only the visible area is illuminated. When zoomed in, the light circle adapts automatically: the higher the magnification, the smaller the light circle. This prevents the possibility of drying or burning exposed tissue, outside of the actual field of view.


Invisible Becomes Visible

Leica M720 OH5 Intraoperative Fluorescence

Fluorescence technologies provide intraoperative information to the surgeon and OR team, directly through the microscope eyepieces or on a monitor. The information gained allows the surgeon to make faster progress in work, increase surgical precision, and improve patient outcomes. Switching between white light and fluorescence mode requires only the push of a button on the hand grip or foot control. The Leica M720 OH5 is well prepared for new and future types of surgical fluorescence, with a selectable third fluorescence mode.

Oncological Fluorescence Leica FL400* intracperative module is


used in conjunction with 5-ALA fluorescent agent to show t umor cells, and thus enables much higher accuracy with tumor resection.



Malignant glioma, blue-violet light mode

Vascular Fluorescence Leica FL 800 [×] intraoperative videcangio-

graphy module is used in conjunction with ICG fluoresc ent agent and allows surgeons to see blood flow through vessels in real time during surgery.



IC G in jection aft er 9 se cond s: veno us view

Surgical fluorescence: The study of fluorescence microscopy has a long tradition at Leica Microsystems, dating back to the beginning of the 20 th century. An indispensable component in biological research, fluorescence science is now an inte gral part of the surgical workflow to improve the patient's quality of life.

*Please check the status of Leica PL400 and Leica PL800 requilatory approval for your country with your local Leica Microsystems represent ative.

Leica HD C100 Camera

The Leica high-definition medical-grade camera delivers bright, sharp pictures and videos, and features an innovative one-touch control button for easy use.

Leica Vide o Adapters

Leica HD video adapters offer intraoperative fine focus and manual or remote control, to always achieve crisp and clear image quality in documentation.

Integrated HD Monitor The Leica M720 OH5 features a built-

in, movable video monitor arm with three rotation axes and an inclination axis to easily maneuver the large videoscreen into the perfect position for all viewers.

HD Documentation Systems The Med X Change HDMD® 1080p or

720p is a user-friendly digital recording system for the surgical environment. The 1080p version records videos in Full HD and detects ICG automatically. Image and video files can be transferred to a USB, external hard drive or wirelessly to an Apple® device within seconds.

3D Documentation System

The TrueVision® 3D Surgical* system combines 3D visualization and guidance software applications focused on improving accuracy, efficiency, and outcomes for both surgeons and patients. AGE ACLUVISIT

SING

 \ast Plea se check the status of Tru eVision $^{\otimes}$ product availability with your local Le ica Microsystems representative .

Leica M720 OH5 Viewing for the OR Team

Leica M720 OH5 OpenArchitecture[™] allows for easy upgrades of rapidly evolving imaging technology. User-friendly operation ensures easy recording and editing of videos and photos for presentations, teaching, or medical records.

IGS Integration

The Leica M720 OH5 includes mechanical and electronic interfaces to accept and easily integrate commonly used imageguided surgery (IGS) systems and their tool tracking cap abilities.

Leica DI C700

The Leica DIC 700 dual in aging color module allows the surgeon to inject data directly into the eyepiece, from external and internal sources, such as MRI, CT, IGS, endoscopes and Leica FL800 video sequences. The data are displayed with the highest resolution and contrast currently available.

Leica M720 OH5 Positioning Flexibility to Suit Surgeon and OR Staff

14

The Leica M720 OH5 provides ultimate positioning flexibility with the highest overhead clearance and longest reach of any surgical microscope on the market. Superior reach and a compact footprint, give the surgeon positioning flexibility to place the microscope wherever it is most beneficial for the surgery. Alternatively, the Leica OHC5 ceiling mount option optimizes performance in space-restricted ORs.

Freedom of Positioning

COLOR MICH

None of Concession, Name

Efficiency in Work: The compact base of the Le ica M7200 H5 creates a smaller footprint, yet provides superior reach and a mple overhead clearance to work in comfort during any surgical case.





Precise Movement





0--



Leica M720 OH5 Superior Maneuverability

Optics Carrier Tilt [1]: The improved inclination angle combined with the most compact optical system provides the surgeon with unmatched comfort and gives much more flexibilityfor transsphenoid al and posterior fossa cases.

Optics Carrier Lateral Movement [2]: With the long est range of lateral movement available today, the surgeon can easily achieve the most challenging lateral approaches.

Erg oLock[™][3]: The main surgeon's bin ocular tub e can be easily locked in five defined position s, ensuring stability of an individual's selected bino cular position, especially when using the mouth switch.

Mouth Switch (Optional) [4]: The ergo nomicallydesigned mouth switch allows the surgeon to easily position the micro scope while leaving both hands free for surgery.

Bra kes: Silen t, high -pr ecision electromag net ic Leica OH technology.

Hand Grip [5]: The ergonomic de sign and sturd y, all-metal construction of the hand grip ensure comfort and sta bility when moving the micro scope.

Foot Control (Optional) [6]: For maximum mobility and for fast, easy adjustments, Le ica Microsystems of fers four models of foot controls: cabled or wireless, 12-funct ion or 16-function. The Leica M720 OH5 offers a greatly expanded range of movement in all dimensions, with intuitive functionality and minimal vibration at all magnification levels.

AutoBalan ce [7]

The hard key "AutoBalance" on the stand saves valuable time. With only two pushes of one button, the system fully balances all six axes quickly and accurately.

Intraoperative Auto Balance [8]

A microscope may need rebalancing during surgery due to changing needs for the surgeon's and assistant's positioning. With one push of the AC/BC button, conveniently located above the optical head, the surgeon can rebalance in seconds, even through the sterile drape.







Leica M720 OH5 **Technical Specifications**

The Leica M720 0H57 0H05 surgical microscopes feature innovative Horizontal Optics Technology for more room to work, a Small Angle Illumination system for better depth perception, and an OpenArchitecture[™] platform to integrate the newest imaging technologies such as Full HD.

ELECTRICAL DATA			
Power connection	1600 VA 50/60 Hz	Compact dimensions	Only 72 m m m inim al height from the m ain surgeon's binocular to the objective, with the microscope in a horizontal position
	100 월 (410 영 / -15 영), 120 월 (410 영 / -15 영), 200 월 420 영 / - 15 영), 200 월 430 영 / - 15 영),		Only 232 mmm inimal l ength from the main surgeon's bi nocular
Cofoty close	220 V (+10 % / -15 %), 240 V (+10 % / -15 %)		to the objective, with the microscope in posterior fossa seated
Safety class	Class I		patient position
FICE M720 MICROS	0185	Surface coating	Cove red wi th a ntimi crob ial coa ting (Ag Protect $^{\scriptscriptstyle M})$
	APO OptiChrome [™] -6:1 zoom, m otoriz ed	O PTICA L DATA	
ELECTRICAL DATA Power connection Safety class LEICA MT20 MICROS Magnification Focus Eyepieces Objective Illumination Main light source Emergency lamp Autolris TM BrightCare Plus SpeedSpot TM Binocular tubes Ernol lock TM	Revolutionary new optical concept with horizontal zoom for maximum compactness of the microscope	Magnilication range	1.5× to 17.0× with 10× eyepiec e
Focus	Motorized via multifocal lens, with manual adjustment	Field of view	12.5 mm to 143 mm with 10×eyep iece
	Widefield exertic coeffor execution and opposite assistant, 12.5× for lateral assistant, dioptric	MICROSCOPE CARRIE	R
	settings±5 with adjustable eye cup	Rotation of optics	540°
Objective	APO OptiChrome [™] multi foœ I le ns	Lateral tilt	50° to left / 50° to right
	200 mm to 500 mm va riable working distance through motorized function. with manual override	Inclination tilt	-3 0°to 4150°
Illumination	Continuo usiva diustable illumination liel diameter with gaussian-shaped light distribution; continuously adjustable brightness at a constant color temperature	XY speed	Zoom-correl ated XY spee d
munimation		Balancing	A, B, C, and Daxes are fully automatic, each can be manually balanced
Main light so urc e	High-performance 400 Watt senon and tamp through fiber optic	Intraoperative	AC/BC button for automatic intraop erative re-balancing of the
Emergencylamp	400-Watt xenon arc-lamp on a sepa rate e lectrical system	balancing	A and C axes, and for re-balancing the B and C axes
Autolris™	Built in automatic, zoom-synchronise di liumination field diameter, with manual override and reset feature	Brakes	One brake for A/B axis, one brake for C axis
		Indicator	LED for fluorescence mode status, LED for video recording status
BrightCare Plus	Safety feature for the working distance-synchronized light control		
SpeedSpot™	Dual laser focusing device for fast, precise microscopepositioning		
Binocular tubes	Binocul artu bestvatured to a blo b utterlivergo nomic height adjustment for optimal body position at the microscope; 115° varia ble angle: 0° to 115° range for main surgeon, -55° to +60° for opposite assistant		
Ergo Lock™	Built in lodging device to hold main surgeone binocularitube fixed in five predefined angles 10°, 38°, 68°, 30°, and 115°		

ACCESSO RIES (OPTIO NA L)

Second observ er	Stereo attachment to beam splitter for second observer
Binocular tube	Variable angle 30° to 150° for second observer
Video adapter	Leica Manual Vide o Adapter (MVA), 55mm, 70mm, 107mm focal length, c-munt, with line toque
	Leica Remote Vide o Ada pter (RVA), 55mm, 70mm, 107mm focal knyth, c-munt, with the toque
	Leica Zoom Video Adapter (ZVA), 3:1 zoom, 35mm to 100mm tecal length, c-mount, with the tecas
	Leica NIR Dual Video Adapter (DVA), 60.5 mm, 79.5 mm focal knoth, c-mount, with time tocae
Autofoc us	The Leica Vide o-An alysis Auto focus gives the surgeon mor e precision and greater control by means of keeping the image crisp and clear.
Image inje ction	Le ica LI U700 high resolution, true color dual maging no dule for correlated and non correlated data display, resolution 1024 × 768 pixels, color de pth 24 bit, gray scale 256, contrast >= 1:300, color temperature 2500° – 9000° K
Asepsis	Sterilizable protective glass cover for the objective, sterilizable components for all drive knobs, commercially available drapes epecifically designed for the Leica M 700)
Laser	Laser mi oro man ipulator available from 3rd party
IGS	
Interface / Compatibili ty	Open architecture for IGS systems

Controls	10-function pistol grips for zoom, focus, all-free release of six brakes. Side button to control three user defined brakes
	motorized lateral tilt and in dination (%) and Leica DLC000 functions. Except for the all-free button, all functions are freely programmable.
	Mouth switch for three brakes (XYZ) (optional)
	12-function foot pedal with controls a tranged longitud inally or transversely, 16-function foot pedal with controls arranged transversely, wired or wireless (optional)
	Hand switch (optional)
Integration of documentation	Prepared for i ntegration of video and digital record ing systems Open architecture
Connectors	Numerous built-in connectors for video, IGS, and control data transfer
	12 Volt DC, 19 Volt DC, and AC connections
Carrier for monitor	700 mm l ong
	Flexible arm with 4 axes for rotation and inclination to carry optional video monitor
Materials	All-solid metal construction
Surface coating	Cove red wi th a ntimi crob ial coa ting (Ag Protect™)
Range cantileve r	Max. 1925 mm
Load	Min. 8.0 kg and max. 11.7 kg of accessories to the microscope
Weight	Approx, 310 kg as a fully conligured system
Storage dimensions	1945 mm (height)× 1395 mm (width) × 830 mm (depth)
AMBIENT CONDITION	IS
In use	+10° C to +40° C (+50°F to +104° F)
	· ,

FLUORESCENCE* (OP TIONAL)

Vascular Huorescence	Uptional Leca H2001 savailable in the USU, BU, and most other countries
Oncological	Uptional Letta Hukoonsavariable in the BU, and some other
Nuorescence	countries
* Please check the status	of regulatory approval for your country with your local
Leica Microsystems repre	esentative.

LE OH5 FLOOR STAND

Туре	Floor stand with six electromagne tic brakes
Base	720 mm \times 720 mm with four 360° ro tatable casters of 130 mm diameter e adh ; one central single step foot brake
Balancing	"No brake re lease" Auto-balance
	One button / two pushes for complete, a utomatic balancing of stand and optics
Intraoperative re-balancing	AC/BC button for automatic intraoperative re-balancing of AC axis and BC axis
Swing arm	Paten ted a dva nœ d movement system for perfect balance in six axes, vibra tion-d issipa ting technology
Floor stand control unit	New generation touch panel technology. The latest electronics control for the continuous operation of all motorized functions and illumination intensity. Data displayed via LCD. Built-in BrightCare Plus technology for working distance synchronized illumination control. ISUS [™] Intelligent Setup System, menu selection based
	on unquesoftware for userspectic conjiguration, with built in
	electronic auto-diagnosis and user support. Software-in dependent hard keys for illumination and autobalancing; indicator for main / badup illumination and fluorecom or mode. Open architecture for future software developments.
Light source	400-Watt dual xeno n arc-l amp illumi natio n system and built-in, automatic (after notice), I amp quick chang er

In use	+10° C to +40° C (+50°F to +104° F)
	30 % to 95 % relative humidi ty
	500 mbar to 10 60 mbar atmospheric pressure
Storage	-40 °C to +70° C (-40° F to +1 58° F)
	10 % to 100 % relative humidity
	500 mbar to 1060 mbar atmospheric pressure

LIMITATIONS OF USE

The Loics M720 OHS surgical microscope may be used only in closed rooms and must be placed on a solid floor. It may not be used in Ophthalmology.

CONFORMITY


CE

Council Directive 33/42/EBC on Modical Devices (MDD) and its anonements. Classification: Class I, in compliance with Annex IX, rule 1 and rule 12 of the directive. Medical Electrical Equipment, Part 1: General Requirements for Safety IEC 60601-1; EN 60601-1; UE0 601-1; CAN/CSA-C22.2 NO. 601.1-M90. Electroma gnetic compatibility IEC 60601-1-2; EN 60601-1-2. The Medical Division, within Leica Microsystems (Schweiz) AG, holds the mana gene nt system catilicates for the international standards ISO 9001, ISO 13485, and ISO 14001 relating to quality management, quality a ssurance and environmental management.


Ap ple is a trademark of Apple Inc., reg ist ered in the U.S. and other countries.
HDMD and Med X Change are trademarks of Med X Change Inc., regist ered in the U.S. and other countries.


METRx is a trademark of Medtronic Inc., registered in the U.S. and oth er countries.
True Vision is a trademark of True Vision 3D Surgical Inc., registered in the U.S. and other countries.

Leica M720 OH5




Leica M720 OHC5







A Paradigm Shift in Vision, Comfort, and Flexibility



www.leica-microsystems.com



The fruitful collaboration "with the user, for the user" has always been the foundation of Leica Microsystems' innovative strength. On this basis, we have developed our five corporate values: Fioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement.

MEDICAL DIVISION

What does a surgeon expect from an outstanding surgical microscope? Sharp, clear images, and a modular system aligned with the surgeon and OR staff needs.

Innovations for your practice

From the first surgical microscope with widefield optics in the 1980s to the first microscopes with Horizontal Optics and with LED illumination, Leica Microsystems has been at the forefront of innovation in the development of surgical microscopes.

HD video, fluorescence and refinal viewing systems also demonstrate the continued innovative nature of the Leica team. We strive to provide the surgeon with leading edge technology to enhance performance, surgeon comfort, and patient outcomes.

Leica Microsystems – an international company with a strong network of worldwide customer services:

Active worldwide	Tel.		Fax	
USA - Buffalo Grove/Illinois	+1	800 248 0123	+1	847 405 0164
Canada - Concord/Ontario	+1	800 248 0123	+1	847 405 0164
Australia - North Rydo/NSW	+61	2 8870 3500	+61	2 9878 1055
Austria Menna	+43	1 486 80 50 0	+43	1 486 80 50 30
Belgium Diegem	+32	2 790 98 50	+32	2 790 98 68
Denmark-Ballerup	+45	4454 0101	+45	4454 0111
France - Nanterre Cedex	+33	811 000 664	+33	1 56 05 23 23
Germany- Wetslar	+49	64 41 29 40 00	+49	64 41 29 41 55
haly. Milan	+39	02 574 861	+39	02 574 03392
Notherlands - Rijswijk	+31	70 4132 100	+31	70 4132 109
Portugal Lisbon	+351	21 388 9112	+351	21 385 4668
Spain Barcelona	+34	900 210 992	+34	93 494 95 40
Sweden - Kista	+46	8 625 45 45	+46	8 625 45 10
Switzenland - Heerbrugg	+41	71 726 34 34	+41	71 726 34 44
United Kingdom - Milton Keynes -	+44	800 298 2344	+44	1908 246 312
China - Hong Kong	+852	2 564 6699	+852	2 564 4163
- Shanghai	+86	21 6387 6606	+86	21 6387 6698
Japan - Toloyo	+81	3 5421 2800	+81	3 5421 2896
Korea - Secul	+82	2 514 65 43	+82	2 514 65 48
Singapore	+65	6779 7823	+65	6773 0628

CE

10 M1 721 0en/00 + Copyright ObyLeica Microsysteme (Schweiz) UG, Medical Elvision, CH 3485 Hoerbrugg, 2013 + Print edin Switzerland – N.2013 – gale da – Subject to modificatione. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH.