

elma clean 65 (EC 65)

NEUTRAL CLEANING CONCENTRATE FOR LABORATORY AND WORKSHOP

Description

elma clean 65 (EC 65) is an universally applicable mild neutral cleaning concentrate with corrosion inhibitor. It is suitable for all cleaning applications which require very gentle treatment of the items to be cleaned. *elma clean 65 (EC 65)* has been developed in particular for the use in an ultrasonic bath. It is suitable for the cleaning of

- laboratory equipment and instruments and
- laboratory work pieces.

elma clean 65 (EC 65) is excellent for work pieces made of sensitive materials and

- glass,
- plastics,
- metals,
- rubber.

It removes contamination such as:

- dust,
- finger prints,
- pigments,
- light greases and oils.

Application and dosage

Use 2 – 5 % in an ultrasonic bath; at 20 – 80° C. After cleaning, rinse thoroughly and dry.

Safety recommendations

elma clean 65 (EC 65) is classified as hazardous according to the regulation (EC) No 1272/2008 [GHS] (serious eye damage). Observe also with respect to this the hints indicated in the Safety Data Sheet and always handle chemicals with care.

Physical-chemical characterisation

Density: 1.04g/ml, pH ~7 in original state. Ingredients according to Annex VII, A, EC-Regulation 648/2004 (detergents): 5-15% anionic surfactants, 15-30% non-ionic surfactants, <5% amino-soaps, <5% phosphates.

Disposal

The surfactants in our product meet the criteria for biodegradation as laid down in Annex III of the Regulation (EC) No 648/2004 on detergents.

The cleaning bath can be fed into the public sewage system; observe the local pH limit values and make sure that the contamination contained complies with the local sewage regulations.

European waste code: 20 01 29*, „detergents containing dangerous substances“.

Volumes, storage and transport

Available volumes: 1 litre bottle; 10 litre and 25 litre HDPE-can.

Store in closed original container only in a cool place at a temperature between -5°C and +35°C, protected from heat and direct solar radiation.

Shelf life: 5 years from date of production (see stamp on label).

Classification for all means of transport: no hazardous material.

Elma Schmidbauer GmbH

Gottlieb-Daimler-Str. 17, D-78224 Singen

e-mail: info@elma-ultrasonic.com

17.02.2015 / P:\Produkte\PROInfo\Produkte\englisch\word_GB\p_EC 65_GB.docx

Tel.: +49 7731 882-0

Fax: +49 7731 882-266

www.elma-ultrasonic.com

