







USB port for updates enables direct access to the LATEST FIRMWARE DEVELOPMENTS

## TECHNICAL DATA

/// C-MAG HS 7 control, Ident. No. 0020002694

Max. stirring quantity	20
Motor rating output	9 W
Speed range	50 – 1500 rpm
Heating output	1000 W
Heating temperature range	RT – 500 °C
Heat control	Turning knob

Max. stirring bar length	30 – 80 mm
Adjustable safety circuit	100 - 650 °C
Set-up plate material	Glass ceramic
Set-up plate dimensions	180 × 180 mm
Dimensions (W $\times$ D $\times$ H)	220 × 354 × 88 mm
Weight	4 kg

## IKA Works (Asia) Sdn Bhd

Lot PT6445, Jalan Industri 3/4, Rawang Integrated Industrial Park, 48000 Rawang, Selangor, Malaysia Phone: +60 3 6099-5666, Fax: +60 3 6092-0193 eMail: sales.lab@ika.my



EN



THE NEW C-MAG HS 7 CONTROL

## **C-MAG HS 7 control**

/// The new magnetic stirrer

The C-MAG HS 7 control magnetic stirrer features a square ceramic heating plate. It can stir up to 20 l ( $H_2O$ ) and has a speed range of 50 to 1500 rpm.

The lifetime warranty, offered with this product, reflects our confidence in the quality of this new magnetic stirrer. Because of its reliability, scientists are now able to concentrate more intensely on their work and experiments.

IKA uses hardened glass for the control panel not only due to its particular safety and resistance to chemicals, but also to introduce an entirely new look to the range of magnetic square top hotplate stirrers.









The magnetic stirrer with a lifetime warranty:



3