

Chlorine Release Tablets

Steadichlor

Steadichlor provides practical solutions for applications in chlorination, disinfection and sterilisation.

Steadichlor tablets offer a simple and effective means of water chlorination without the need for complex dosing equipment. Steadichlor tablets dissolve in water to give a slow steady chlorine release. The dissolving time is approximately four hours in flowing water or 12-24 hours in static water. Steadichlor tablets are prepared from calcium hypochlorite. In use, the tablets do not contribute any substances not normally present in treated mains water.



Chlorination of Water

Steadichlor tablets offer an ideal means of chlorinating water in tanks, wells or small reservoirs, and waste water in tanks or effluent streams. The tablets should be used in accordance with standard chlorination practice depending on the water to be treated and the purpose of treatment.

Where chlorination is to be carried out to disinfect the water, a free chlorination residual of at least 1 mg/l should be achieved in the treated water. Note that it will normally be necessary to add a higher chlorine dose to satisfy the chlorine demand of the water.

Volume of water to be treated		No of Tablets to give 1 mg/l chlorine dose Steadichlor - NS
Litres	Gallons	
14000	3700	1
28000	7400	2
56000	14800	4

Control the dosage of Steadichlor NS Tablets
Find out about our Chlorine Test Kits

www.palintest.com



Disinfection of Water Mains and Tanks

The Steadichlor system is established as one of the simplest means of disinfecting new or repaired water mains and tanks. A feature of the Steadichlor system is the unique Palintest Chlorocol test, especially developed for the control of mains disinfection procedures.

Pipe Diameter		Spacing between each tablet (metres) required to give minimum 30 mg/l chlorine dose
mm	in	
80	3	-
100	4	50m
150	6	25m
200	8	13m
250	10	8m
300	12	6m
400	15	3m
450	18	2m

Palintest Test Kits

Palintest offer a range of test kits for testing chlorine

ChloroSense - Unique Chlorine Measurement

Unique to Palintest, ChloroSense uses a unique and patented electrochemical sensor to measure free and total chlorine simultaneously in a potable water sample.

Chlorometer/ Chlorometer Duo

Measures the level of chlorine disinfectant in drinking water distribution (0 - 5 mg/l) and other low-level chlorine residual applications. The Chlorometer Duo also offers high range testing (0 - 250 mg/l).

We also offer a simple pocket kits to perform the Chlorocol Test in a variety of applications. You can find more information on our test kits by visiting our website (www.palintest.com).



Ordering Information

Product Code	Pack Size
WT 122, Steadichlor 'NS'	2.5 Kg (approx 125 tablets)
WT 123, Steadichlor 'NS'	5 Kg (approx 250 tablets)

Palintest Ltd
Palintest House, Kingsway, Team Valley,
Gateshead, Tyne & Wear, NE11 0NS,
England
Tel: +44 (0) 191 491 0808
Fax: +44 (0) 191 482 5372
Email: sales@palintest.com
Website: www.palintest.com

Palintest USA
21 Kenton Lands Road, Erlanger,
Kentucky 41018, USA
Tel: (859) 341 7423 **Toll Free:** 800 835
9629
Fax: (859) 341 2106
Email: info@palintestusa.com
Website: www.palintestusa.com

Palintest Australia/Asia Pacific
1/53 Lorraine Street,
Peakhurst Business Centre,
Peakhurst NSW 2210, Australia
Tel: +61 - 1300 131516
Fax: +61 - 1300 131986
Email: palintest@palintest.com.au
Website: www.palintest.com.au

Palintest China
Room 1601, Kun Tai International
Mansion, 12B Chao Yang District,
Beijing, 100020, PRC
Tel: +86-10 5126 1868 - 809
Fax: +86-10 58790155
E.mail: fred.fan@palintest.com
Website: www.palintest.cn

A HALMA COMPANY

